

CaST II

User's Guide

Version 5.0.50

Prepared by
Information Management Services, Inc.

In cooperation with & supported by
The Centers for Disease Control and Prevention

February 2003

Table of Contents

OVERVIEW	1
Introduction	1
What is CaST?	1
What can CaST do?	1
Where does CaST come from?	1
User's Guide Layout	1
General Layout	1
Technical Specifications	2
Technical Assistance	2
System Overview	2
Changes for Version 5.0.50	2
Network Capable	3
User-defined Fields	4
Query Features	4
MDE Export Feature	4
INSTALLATION	4
Hardware/Software Requirements	4
Running the Installation	5
Installation Types Supported	12
Workstation Only -- Standalone	12
Run from Network -- Workstation Program on Server	12
Network Data -- Workstation with Networked Data	13
Administrator Responsibilities	13
New Programs & New CaST Installations	13
Administrative Setup	14
Upgrading the Databases for CaST 5.0.50	16
STARTING CAST	18
User Name and Password	18
Logging On	18
Logging Off	18
Changing Your Password	19
SYSTEM LAYOUT	19
Main Menu	19
Menu Bar Items	20
File Menu	20
Search Menu	21
Report Menu	21
Tools Menu	22
Administration Menu	22
Help Menu	23
DATA ENTRY	23
Screening Cycles Revisited	23
Cycles and Procedures	23
Required Fields	24
Short Cuts	24

Data Entry Field Types	25
Moving from Field to Field	26
Saving Data	26
Tab-delimited Screens	26
Entering New Patient Data	27
Baseline Information Data Entry	28
Baseline Information	28
Contact Information.....	29
Cycle/Procedure Data Entry	30
Adding Cycles and Procedures.....	30
Cycle Patient Information.....	30
Cycle Procedures.....	32
Cycle Treatment Information.....	33
FIND PATIENT/PATIENT STATUS.....	34
Find Patient Dialog.....	34
Patient Status.....	35
Patient Status.....	35
Working with the Patient Status Screen	36
QUERIES.....	37
Introduction to Queries	37
Retrieving a Saved Query	38
Working with Query Results	39
Creating a New Query.....	42
Creating a New Query -- Getting Started.....	42
Selection Criteria.....	42
Fields Selection.....	45
Sort Options	46
SQL Tab.....	47
Modifying an Existing Query	47
Advanced SQL.....	48
STANDARD REPORTS	48
What are Standard Reports?.....	48
Running a Standard Report.....	49
Modifying an Existing Standard Report.....	49
PATIENT REMINDERS	50
What are Patient Reminders?	50
How to Run Patient Reminders	50
Working with Patient Reminder Results.....	51
How Patient Reminders are Created	53
OTHER REPORTS	54
Overdue Reports.....	54
What is an Overdue Report?	54
How to Run the Overdue Report.....	55
Duplicate Patient Report	55
Finding Duplicates in the Database	55
Running the Duplicate Patient Report	56

MDE EXPORT	56
Creating the MDE Export File	56
Running the MDE Edit Check Program	57
Compressing the MDE File for Submission	59
WORKING WITH DECENTRALIZED DATA ENTRY	60
Import/Export Functions	60
What are the Import/Export Functions?	60
Export Data to the Central Site	60
Import Data from a Data Entry Site	61
Restoring/Creating New Site Databases	64
Creating a Restored or New Site Database	64
Restoring a Site's Database	64
Creating a New Database for a Site	65
Exporting/Importing Provider Sites	66
Using the Provider List Export and Import Functions	66
Exporting Provider Sites	66
Importing Provider Sites	67
ADMINISTRATION	68
User Setup	68
Provider/Enrollment List Setup	70
Data Entry Sites	72
Layout Files for the Procedure Data Entry Screens	73
Database Password/Encrypt Editor	75
Format Library	76
What is a Format Library?	76
The User-defined Format Library	76
Using the Format Editor	77
User-defined Fields	79
Introduction to User-defined Fields	79
User-defined Fields Setup	79
BACKUP/RESTORE/COMPACT	81
Backing up the Database	81
Restoring a Backup of the Database	82
Compacting the Database	83
FREQUENTLY ASKED QUESTIONS (FAQ)	84
Known Issues and Comments from current Users	84
Is there a billing component in CaST?	84
How can I contact other CaST Users?	85
What do 'Old CBE - Normal' and 'Old CBE - Abnormal' mean?	85
What are the 'Old Location' fields?	85
Why doesn't CaST support the use of numbers for data entry?	86
Why 'CaST II' for the directory and logon screen?	86
How can I print the MDE summary and error files?	86
Why doesn't the cycle number database field match the record location in CaST?	87
How do I delete a procedure?	87
How can I get an electronic list of my providers for use in other applications?	87

TECHNICAL SPECIFICATIONS 88

Database Structures..... 88

Database Conversion..... 88

MDE Export Specifications 88

APPENDIX I - III

Overview

Introduction

What is CaST?

CaST is an acronym that stands for **C**Ancer **S**creening and **T**racking System, and it is designed for all levels of users, from novice to expert. CaST is compatible with Windows 95, 98, NT, Me, 2000 and XP and uses standard Windows layout for its screens and has menus for access to its features. If you are not familiar with basic Windows Operations, please refer to your Windows User's Guide *before* using CaST.

What can CaST do?

CaST allows the user to track women through the length of the study and collect information on screening and diagnostic procedures done for breast and cervical cancer. To assist in the tracking of women, CaST has the ability to generate reports, queries, and client/patient reminders. Secondly, it provides a means to report the Minimum Data Elements (MDEs) to the Centers for Disease Control and Prevention (CDC).

Where does CaST come from?

CaST has been provided by the CDC through contract with Information Management Services, Inc. (IMS) for Programs in the National Breast and Cervical Cancer Early Detection Program (NBCCEDP) to offer a user friendly environment to collect data.

User's Guide Layout

General Layout

The CaST User's Guide is written in a style that makes it easy for you to find the information that you want. Each chapter of this guide corresponds to a major component in CaST. For example, there are separate chapters for Data Entry, Queries, Reports, and Administration. Within each chapter, different sections explain how to use each major component of

CaST. For example, the Queries chapter contains a section on "Working with query results".

The User's Guide can also be accessed within the system by selecting 'Contents and Index' from the 'Help' menu.

Technical Specifications

At the end of the printed copy of the User's Guide, there is an appendix which contains technical specifications, including data structures.

Technical Assistance

To obtain technical assistance on CaST, please contact the CaST Technical Support Team at Information Management Services, Inc. You can call the team at (301) 680-9770 or email them at casthelp@nbccedp.org.

System Overview

Changes for Version 5.0.50

CaST, Version 5.0.50 is an upgrade from version 4.1.1, which was distributed in April 2001. Existing users of version 4.1.1 must also run an update 'script' to convert the CaST databases. Please consult the Upgrading the Databases for CaST 5.0.50 topic in the "Installation" section of this document.

The major changes from version 4.1.1 to 5.0.50 are as follows:

- The system and database were modified to implement the new MDE 5.0 data definition. The following changes were made to incorporate the 5.0 changes:
 - Four race fields were added to the baseline table and data entry forms.
 - The existing 'Pap Smear' procedure was renamed to 'Pap Test 1991'. A new procedure was added that is called 'Pap Test', that will collect the Pap test results using the new 2001 Bethesda Reporting system.
 - A new field called 'SpecimenType' was added to the database, along with the cervical procedures data entry

area, to collect Specimen Type for the new Pap Test procedure.

- A new cervical procedure, "HPV Test", was added. This procedure name is a shortened name for the "Digene hc2 HPV Test". This procedure will currently be ignored in the MDE export.
- The baseline data entry form now has the enrollment information and comments section at the top of the form. The chart number and CD number fields were also moved to the top, but are not part of the tab order in data entry.
- The system will now allow for more fields to appear on one screen. Each user can maximize the screen (or just resize) to fit more fields on the viewing screen. The system will remember each user's preference the next time the form is brought up.
- The fields which get exported to the MDE file are now color-coded. The field labels appear in 'red' and the fields are 'pink'.
- The provider list editor now allows for 'inactive' providers. If you uncheck a provider (either enrollment, primary, provider check boxes) and records exist for that type, the system will mark that provider as 'inactive'. 'Inactive' providers will still show up in the data entry screens, but they will appear at the bottom of the list. The 'active' providers will show up first and be used for the 'type ahead' feature.
- Administrator can now password protect or encrypt the database.
- The system no longer uses ODBC to connect to the Access databases. As a result, the Access database can be converted to Access 2000 if needed.
- Crystal Reports is no longer used to generate reports. Users now have the option to print the report directly to a PDF file.
- The functions to handle the transfer of data between decentralized data entry sites and the central site were enhanced.

For a complete list of changes for version 5.0.50, please refer to the document located in the EXTRAS directory on the installation CD-ROM.

The following sections will highlight some of the features in CaST.

Network Capable

CaST is a multi-user system, capable of residing on a network. Although the entire system, including executables could be run off of a network, the preferred method would be having only the data on the network.

User-defined Fields

CaST allows each Program to setup and collect additional fields. There is a user defined table, each able to contain up to 200 fields, for each main table in the database, including the baseline, cycle and procedure tables. These fields will be generated and maintained by the system administrator, and will be available for data entry and for use in all queries in the system. User defined formats can also be created in conjunction with these fields, so that end users will enter data using 'drop down' lists, instead of free text.

Query Features

CaST has query capabilities, which allows the user to execute pre-defined searches, or create searches from scratch. Any field in the database, including user defined fields, can be used in a query. Pre-defined queries include logic based on the MDE feedback audits distributed with the standard MDE reports.

MDE Export Feature

CaST has a built in MDE export function, which allows the user to create a MDE file with a single button click. The MDE export function is automatically updated when changes are made to the MDE definition. The system also allows the user to run the MDE edit check program directly from CaST.

Installation

Hardware/Software Requirements

The minimum configuration recommended for use with CaST is a Pentium 166 IBM compatible microcomputer with 32 Megabytes of RAM. 200 Megabytes free space on the hard drive is recommended, but this could be significantly less or more, depending on the size of the database and if the database is stored on a network.


The operating system needed is Windows 95, 98, 2000, Me, XP or

Windows NT 4.0 with Service Pack 5 or later, with an 800 x 600 or higher screen resolution. To install CaST on Windows 95, 98 (not including Windows 98 Second Edition), or NT, the Y2K Update 2 (for Windows 95 and 98) and Internet Explorer 4.01 with Service Pack 2 or later are also required.


Running the Installation

Your system administrator may have already set up the installation parameters for you. If this is the case, you will be running Setup.exe from a network location. Simply verify that the parameters are correct, and contact the administrator if you think there is an error. If you are specifying the parameters yourself, follow these instructions.


To install CaST 5.0.50, first close any open programs or utilities. Place the installation CD labeled "CaST 5.0.50" into your CD-ROM drive and run the "Setup.exe" program located in the "Install" directory. The first screen is a welcome screen; click 'Next' to continue.




You will be prompted to choose a type of installation: Workstation only, Run from Network, or Network Data. Highlight your choice and click 'Next'.




The next screen asks for a destination directory for the program. The default is "C:\Program Files\Cast II". It is recommended, but not required, that the files be installed in the default directory. To select a different directory for CaST, click 'Browse' and locate or create the new directory. Click 'Next' to continue after verifying the installation directory.




Next you will be prompted to supply a directory for the CaST support folders and files. The default location is "C:\Cast Files\". You can change this location as well by clicking 'Browse', but it is generally recommended that you choose a location on your local hard drive. Click 'Next' when ready.




On the next screen you must provide the locations for the Data, Format, and Layout files. The default location for all is a subfolder called "Data" in the support files folder you set up previously. (For Network Data installations, this will be a location on the network.) It is recommended that you use the same location for all files, and for this reason, using the 'Browse' button next to the Data location field will change the other locations as well. However, using the other 'Browse' buttons gives you the opportunity to specify separate folders for the three types of files. When you have selected all three locations, click 'Next'.



The next screen will prompt for the folder in which to put the program icons. If you are using Windows 95/98 or NT 4.0, this will be a folder on your Start menu. If you are using Windows NT 3.51, a program group of this name will be created. "CaST II" is the default choice. If you want to change it, highlight an existing folder from the list, or type in the name of a new one. Click 'Next' when ready.



If you are performing a Workstation only installation, the installation program will now search your hard drive for existing data files from previous versions of CaST. Click 'Next' to begin the process. The program will look for "user.fmt", "castmain.mdb", and "users.mdb" in your existing data directory, if any, and will continue to search your entire hard drive until the files are found. When it finds one of the files, it will prompt you to decide whether to copy the old file to your new data location. Verify that the correct file is being used before clicking 'Yes'.



Before beginning to create new directories and files on your computer, the installation program will display the information you entered on the previous screens for verification. If any information looks incorrect, click 'Back' until you reach the screen where you can change it. If you are satisfied with the installation parameters, click 'Next'.

A progress meter will track the installation process.

If you are a new program installing CaST for the first time, you will be prompted to enter your **program** and **site codes**. This is a very important step; it determines how each patient within the database is numbered.

program code – The two-letter abbreviation for your program.

site code – This code is used for decentralized sites. If you are not a decentralized site, enter 001. For a decentralized site, put in the code you obtained from your central site.

Click 'Finish' when the final screen of the installation program is displayed. You are now ready to use CaST.

Important Note: System administrators, please read the following topics concerning administrative installation instructions and upgrading from previous versions of CaST:

- New Programs & New CaST Installations
- Administrative Setup
- Upgrading the databases for CaST 5.0.50

Installation Types Supported

Workstation Only -- Standalone

A standalone installation is one in which a single machine will be running CaST and storing all data on the local hard drive. Everything is contained on the single machine and thus there is a single point of administration. Regular backups of the data must be made to tape or other media to ensure data availability in case of a hard disk failure. This type of installation is suited for small programs that have only one person in charge of the data entry and all reporting, or that do not have a network. It is the least flexible of the installation types.

Run from Network -- Workstation Program on Server

A network installation refers to having all of the programs and data residing on a network server or a machine acting as a network server. In this installation all of the programs and the directory structure are contained under a single program directory created and managed on the server. This installation is the preferred type of installation from an administrative point of view. All of the files are located in a central directory and there is only one copy. Only one set of files need to be replaced when upgrading. Individual users will run the executable from the shared network point, and will therefore always be running the current version. The drawback to this type of installation is reduced speed. Because of the limitations of most networks, the program will execute more slowly when run from the network, as opposed to from the hard disk on a local machine.

Network Data -- Workstation with Networked Data

A workstation installation with networked data is a compromise between a standalone and a fully networked installation. In this installation, the program and all of its support files are stored on the local machine, while the database is stored on the network server. This type of installation should provide faster program execution than the full network installation. All files except the data are stored locally; therefore, less data is transferred while the program is running. This type of setup allows multiple workstations to access the data at the same time, and also allows the server to handle the backup of the most important files, the data. A drawback to this installation is the additional work necessary for the administration of the program. Since every machine has a copy of the program and associated support files, each machine must be upgraded when a change is made to either the program or a support file.

Administrator Responsibilities

New Programs & New CaST Installations

As an administrator, there are important steps you must follow in the installation and setup of CaST. You must first determine which installation type your program will utilize. Once this is determined, you must make sure that the first installation is performed according to the directions in the following section. Finally, you must configure the system options, including user setup, provider/enrollment list setup, user-defined fields, format library, and 'layout' files.

When installing CaST for the first time as a new Program, if a networked type of installation is being used, there are differences between setting up the workstations and the initial installation. (For Workstation Only installations, there are no further instructions.) Follow these steps before conducting the workstation installations in order to avoid conflicts.

For "Networked Data – workstation with networked data" installation, you must first install CaST to a local machine. First, run a Workstation Only installation, and use the default location – on the local hard drive. Next, copy **all files** from the Data directory up to the location on the network. You will now be able to run the installation on other workstations and choose "Networked Data" as the installation type. The installation program will prompt you for the location of the data directory. Use the UNC (i.e. \\server1\nbcc\cast\data) path name for this location. To use CaST on the original machine, you need to uninstall and then reinstall it, using the

"Networked data" option.

For "Run from Network – workstation program on server" installations, you must first make sure all program files and databases are located on the network. You could install CaST directly to the network, but this is not recommended, because a subsequent uninstall would remove all the files from the network. A 'safer' alternative would be to first run a Workstation Only installation, and use the default location – on the local hard drive. Then copy the entire "CaST II" directory, along with the Data directory, up to the network. Run the installation on other workstations and choose "Run from Network" as the installation type. The installation program will prompt you for the location of the data directory and the location of the executable. Use the UNC (i.e. \\server1\nbcc\cast\) name for both of these locations. You will need to uninstall and then reinstall CaST on the original machine, using the "Run from Network" option.

The program is installed with only one active account, **Admin**, with no password. This account has all the ability to create users and change any of the program settings. **It is highly recommended that you first set a password on this account.** The password can be set by selecting the "File | Change Password" menu item and entering a new password in the dialog box.


Administrative Setup

System administrators can set up a simplified installation process for their users, in which the desired choices and file locations will be filled in by default.

This Administrative Setup should also be run by those programs performing Network Data or Run from Network installation. Administrative Setup will update the format and layout files necessary to run CaST 5.0.

To begin, run CaST's setup.exe program, including the **/a** switch on the command line. For example, if D is the letter of a CD-ROM drive with the CaST installation CD inserted:


1. Click the 'Start' button at the lower left of your screen.
2. Click 'Run'.
3. Type "D:\Install\Setup.exe /a".



The setup program will start in Administrative Install mode. Click 'Next'.

You will be prompted to select a directory to which the administrative installation files will be written. Most likely, this will be a location on your network which is accessible by all users. Type in the field or use the 'Browse' button to set the directory.

Follow the installation process as described in Running the installation, setting the parameters that every user should have.



After the installation information is entered, the next screen allows you to determine what information will be written to an INI file that will accompany setup.exe. If you do not check the option to create the INI file itself, the information you have entered will not be available to users when they run

the installation program. For each of the options that you check, the appropriate information will be included in the INI file, and used as the default parameters during installation.

Individual users may still change the parameters as necessary, unless you check the 'Quick Install' box. Quick installation uses all the parameters you have specified without giving the user the opportunity to edit them. Users are presented with the verification screen, and may choose either to accept the parameters, or cancel the installation.

The final screen prompts you to specify which program and data files you want written to the directory you chose at the beginning. These files will be written immediately. If you chose the "Run From Network" installation, these files include updated format files, layout files, and all program files.


The installation files, including the INI file, will then be copied to the network location you specified on the first screen. Instruct users to execute the setup.exe program located in this network directory in order to take advantage of the easy installation process.

Upgrading the Databases for CaST 5.0.50

Once you have upgraded the workstations, you will need to run a conversion utility to upgrade the databases. This will only need to be run once, and should be performed by the system administrator.

There are two files which you will access on the CD-ROM. These files are located in the "Utilities" directory and are called "CaSTUpdate.exe" and "CaSTUpdate40.cii". (Note: if you previously downloaded and ran the CaST 4.11c patch, use the "CaSTUpdate401.cii" file located in the "411cpatch" subfolder.)

To run the update, make sure there are no users in the system (for networked data users). Using Windows Explorer, double click on the "CaSTUpdate.exe" file. The following screen will appear:



To perform the database update, do the following:

- Click on the "Browse" button.
- A standard "Open file" dialog box will appear, choose the "CaSTUpdate40.cii" or "CaSTUpdate401.cii" file (see note above), and click 'Open'.
- The file name will now appear in the "Input Script" box. Click on the "Update" button to start the database upgrade.
- This update utility will first back up the two databases. It will choose the default "TEMP" directory for your machine (i.e. c:\windows\temp). Click "Ok" if appropriate or change the location for the backups. Note: these files may already exist on your machine (from previous database upgrades) – you can overwrite these files.

The script will then display information in the "Messages" area. The script should take no more than 10 minutes, and may only take a few seconds. The time depends on the size of your database and the speed of your machine. When completed, a message "Update process completed" will appear. Click "Ok" and "Close" the update window.

The databases are now upgraded, and you can run CaST 5.0.50.


Starting CaST

User Name and Password

To use CaST, you must be assigned a user name and password. CaST reads a special database called the Administrator database, which contains a list of the current users and their privileges on the system. The user names for CaST are assigned by your system administrator.

Logging On

The first thing that CaST does when it starts is display the logon box and prompt you for your user name and password.

A screenshot of a 'Database Login' dialog box. The dialog has a title bar with 'Database Login' and a close button. Inside, there is a label 'Database:' followed by the text 'CaST II'. Below this, there are two input fields: 'User Name:' and 'Password:'. At the bottom, there are two buttons: 'OK' and 'Cancel'.

Type in your assigned user name and your password, then click the OK button or hit the <Enter> key. CaST then verifies that you are a valid user and opens the program. Note: passwords **are** case sensitive.


If you have trouble logging on or forget your user name or password, please contact your system administrator.

Logging Off

When you are finished using CaST, select Exit from the File menu or click the 'x', or close box, in the upper right hand corner of the main menu. CaST then logs you off of the system, exits, and returns you to Windows.

Changing Your Password

To change your password, select Change Password from the File menu on the main menu. The following screen will appear.




A dialog box titled "Change User Password" with a close button (X) in the top right corner. It contains four text input fields: "User Name:" with the text "admin", "Old Password:", "New Password:", and "Verify Password:". At the bottom right are two buttons: "Apply" and "Close".

If you forget your password, please see your system administrator. If you are the system administrator and forget your password, please call the CaST Technical Support Team for assistance. Note: passwords **are** case sensitive.

System Layout

Main Menu

After you enter your user name and password, the following main menu appears:



A screenshot of the "CaST Ver. 5.0" main menu. The title bar shows "CaST Ver. 5.0" and standard window controls. The menu bar includes "File", "Search", "Report", "Tools", "Administration", and "Help". The main area contains five buttons with descriptions:

New Patient	... Enter data for a new patient
Find Patient	... Find a patient within the current database
Standard Reports	... Generate a standard report
Queries	... Generate a user defined listing of the data
Generate MDE	... Create MDE file

At the bottom, a blue banner reads: *National Breast and Cervical Cancer Early Detection Program*

Each section of the system can be accessed using the main menu, either from the buttons or the menu bar. The most frequently used options in the system are available as buttons on the main menu. The options can also be accessed from the menu bar.


You can select a button or menu bar item by either using the mouse or by using a keystroke sequence. For example, to select "New Patient", press the "ALT" and "N" keys at the same time. All buttons in CaST have one character underlined and can be accessed in the same manner. Similarly, to access the menu bar, press the "ALT" key and the corresponding first letter of the menu item. For example, to select the "File" menu, press "ALT" and "F". Once the menu bar is selected, you can use the up and down arrow keys or the underlined letter of the menu item to access the various items on the menu bar.

Note for Windows 2000 users: the underline described above MAY not appear until after you hit the "ALT" key. You can control how these 'hot keys' function by going to the Control Panel, selecting 'Display', then selecting the 'Effects' tab. There is a check box labeled 'Hide keyboard navigation indicators until I use the Alt key'. If this option is checked, the menus in CaST will not have the underline as described above.

Each option on the menu bar is described in the following section, Menu Bar Items.

Menu Bar Items


File Menu



New Patient – identical to the "New Patient" button on the main menu.
Change Password – changes password. See Changing your password.
Provider List – Edit the Provider/Enrollment List.
Export/Import provider lists.
Re-sync – see Provider/Enrollment List Setup.


Page Setup – not activated
Printer Setup – not activated
Exit – closes all databases and exits CaST.

Search Menu




Find Patient by Name – identical to Find Patient button on main menu.
Find Patient by ID, SSN, Disease #, DOB - Find Patient using different fields.
Find Duplicate Patients – runs the Duplicate patient report.

Report Menu



Standard Reports – identical to the Standard Reports button on main menu.
User Defined Query – identical to the Queries button on main menu.
Overdue Report – runs the Overdue Report, which looks for missing data.
Patient Reminders – creates a list and mailing labels for women due for services.

Tools Menu



Generate MDE – identical to the Generate MDE button on main menu.

Backup Database – allows you to Backup the databases.

Restore Database – allows you to restore a previously backed up database.

Compact Database – allows you to compact the database (for efficiency).

Convert Patient Ids – used by decentralized sites for Patients transferring from another site.


Export Data to Central Site – used by decentralized data sites to transfer/send data to central site.

Import Data at Central Site – for programs using decentralized data entry

Format Library – viewing/editing formats used in the system, including those used in user-defined fields.

Administration Menu

For those users with Administrator privileges, the following menu items will be visible.




User Information – user setup.

Program Defined Fields – modifies user-defined fields included in the system.

Site Information – used to track decentralized data entry sites.
Site/Program Editor – used only to change program or site information.
Database Password Editor – used to set the database password or encrypt the database.
Layout file location Editor – sets the file location for the layout files for the procedures data entry grids.
Reset Patient Counter – resets counter used to create Patient ID.
Decentralized Setup – used to create or restore a decentralized site database.

Help Menu



Contents and Index – accesses the online help
Technical Support – displays contact information for technical support
About – displays information about the CaST version, database location, and other useful information

Data Entry

Screening Cycles Revisited

To perform data entry in CaST, you must have an understanding of the definition of the MDE screening cycles. Screening cycles as defined by the MDEs are different than a screening cycle defined in a clinical environment. The user is required to enter data into CaST using the MDE screening cycle definition. If you need help understanding screening cycles, refer to the MDE Data User's Manual or contact your Program's Data Manager.

Cycles and Procedures

Data entry in CaST is comprised of baseline data, along with cycle and procedure information for both breast and cervical cancer screening and follow-up. Each patient in the baseline table has corresponding records in the cycle table, and each cycle has corresponding records in the

procedure tables. Cycle information includes location, referral information, and final diagnosis and treatment results. Most cycles begin with a screening procedure, either a Mammogram, a CBE, a Pap test, or a pelvic exam. If further diagnostic tests are necessary, additional procedures will be entered for that cycle.

Required Fields

Due to the various ways in which Programs enter and perform quality control on their data, there are few edit checks performed during data entry. This includes required fields. In order to add a new patient, the Social Security number or Last Name must be entered. Besides this, and some basic edit checking, the majority of the edit checks will be performed using the MDE Edit Check, the Standard Audit queries, or the Overdue Report.

There are, though, required fields for the MDEs. Fields that are exported to the MDEs are shaded in pink, with red field labels. The MDE Edit Check will also help you determine those fields that were not completed properly. Appendix II also provides the database and data entry fields that are used in the MDE file.

Short Cuts

Throughout this manual you will find helpful 'short cut' keys that will help you navigate the system. Short cut keys are provided to eliminate the need to use the mouse. Although the following can be found in other sections in the manual, they are provided here for quick reference.

- To select a button: press the 'ALT' key and corresponding underlined character. For example, on the main menu, to select "New Patient", press the 'ALT' key and "N" and the same time.
- To Move Forward: press the ENTER or TAB key.
- To Move Backwards: press the SHIFT and TAB keys at the same time.
- To display the values in a drop down field, or list box, press F4. You can then scroll down to the correct value using the up/down keys, or select the value by typing the first few characters (type-ahead feature). Note: you can use the type-ahead feature without

the F4 key.

- To cancel changes: press the ESC key. There is also a 'Cancel' button on most pages. Note, pressing ESC key will only work on certain fields. Also, pressing the ESC key once will cancel the changes made for the current field, pressing it again will cancel the changes made for the entire record.
- To move to the next TAB: press the CTRL and TAB keys at the same time. (i.e., to move from the Baseline Information to the Contact Information tab) Note: if you are in the last field on a tab, pressing the ENTER or TAB keys will also move you to the next tab.
- Moving the cursor from the procedures 'grid' to Workup Planned: to exit the procedures screen, press the F2 Function key to move to the 'Workup Planned' field. F2 will also move the cursor from the 'grid' to the bottom buttons on the 'Find Patient' screen.
- To automatically run a query or report, double click on the item. This will run the query or report as if you selected the item and clicked 'Go'.

Data Entry Field Types

There are several different types of fields in CaST. Examples of each type and a brief description is provided below.

A screenshot of a text input field for a name. The label "Name (Last, First):" is in a grey box. To its right are two adjacent white text boxes separated by a comma, each with a thin grey border.

Text box: used for free text or numbers, only for set length fields

A screenshot of two list box controls. The first is labeled "Ethnicity:" in red text, followed by a white box with a downward arrow. The second is labeled "Race(s):" in red text, followed by two adjacent white boxes, each with a downward arrow.

List box: used for formatted fields, user can only select items in list (press F4 to display values)

A screenshot of a date input field. The label "Date of Birth:" is in red text. To its right is a white box with a thin grey border, containing three underscores separated by slashes, representing the mm/dd/yyyy format.

Date fields: user must enter date using mm/dd/yyyy format. Note: dates

located on the baseline and cycle pages will have the '/' visible before data entry, dates located on the procedures 'grid' will display the '/' as the user types in the date.



Memo fields: used for free text fields, no maximum on number of characters.

Moving from Field to Field

Once information is entered for a field, you can move to the next field by either pressing the Enter or Tab key. To move in the opposite direction, press the "Shift" and "Tab" keys at the same time. You can also move to another field by using the mouse.

Saving Data

Each data entry page has a Save and a Cancel button in the bottom right hand corner of the screen. To post changes to the database, select the "Save" button. This will validate the data and save the changes to the respective tables in the database. For data entry screens with multiple pages, or tabs, of information, the Save button will post changes from all tabs.

If you are adding or modifying data and don't want to save the changes, select the Cancel button. This will ignore any changes to the record and restore the information to what it was before the last Save was performed.

Note: you can also press the 'Esc' key to cancel changes in a particular field or record.

Tab-delimited Screens

All of the data entry screens in CaST implement "tab-delimited" functionality. These screens allow multiple pages of information to be contained in one screen, allowing navigation of the various pages via tabs. Each page on the screen is labeled and can be accessed by either using the mouse or by pressing the 'Ctrl' and 'Tab' keys at the same time.

Each tab has a default size with a set number of fields assigned to each tab. Each user can resize (and even maximize) certain screens to fit more fields in each tab. The system will remember each user's preference for screen size. The example screens provided in this manual are shown at their default sizes, based on CaST version 4.11 and earlier.

Note: During data entry, pressing the 'Tab' or 'Enter' key from the last field on the page will also move you to the next 'Tab'.

Entering New Patient Data

When you select "New Patient" from the main menu, a unique patient ID number is created, never to be used again. This number is created using the program code, the data entry site, and a sequential number. This patient ID, along with the patient name is displayed on the top of the data entry windows. The baseline information screen contains two sections, Baseline Information and Contact Information.

Baseline Information Data Entry

Baseline Information

The screenshot shows a software window titled "Baseline Information Data Entry - - 039001000000148". It has two tabs: "Baseline Information" (selected) and "Contact Information". The "Baseline Information" tab is divided into two sections: "Enrollment" and "Identification".

Enrollment Section:

- Location:** A dropdown menu.
- Date:** A date input field with slashes.
- Chart #:** A text input field.
- CD#:** A text input field.
- Comments:** A large text area.

Identification Section:

- SS Number:** A text input field.
- Alt Id / Type:** Two text input fields.
- Name (Last, First):** Two text input fields.
- Middle:** A text input field.
- Maiden:** A text input field.
- Date of Birth:** A date input field with slashes.
- Ethnicity:** A dropdown menu.
- Race(s):** Six dropdown menus arranged in two rows of three.
- POB City, State:** A text input field followed by a state dropdown menu.
- Country of Birth:** A text input field.

Buttons: At the bottom, there are buttons for "New", "Delete", "Breast", "Cervical", "Save", and "Close".

The Baseline Information section contains baseline information for the patient.

The 'Save' button will post changes to the database. The 'Cancel' button will ignore any changes made since the last Save and restore the record to its previous state. You can use the buttons on the bottom to enter Breast or Cervical information.

You must enter at least the Social Security Number or last name of the patient to save the baseline record and continue to enter cycle information. As you enter information into the Baseline, the system will automatically check to see if this new patient is a duplicate.

Helpful Hints:

- You will notice six 'Race' fields. Only one is required to be completed, but for those women that indicate more than one race, the additional fields can be used.
- "POB" stands for "Place of Birth".

- The drop down category for "Enrollment Location" is populated with those entries in the "Provider List" that have the enrollment box checked.

If there is more than one 'tab' of data entry, click on the next tab or press "Ctrl" and then "Tab" to proceed to view/enter the remaining Contact Information. Note: if you press the 'Tab' or 'Enter' key while in the last field on a 'tab', the system will move over to the next tab. Alternately, you can click and drag the lower right corner of the window to expand it. If the window is large enough to accommodate additional sections from the Contact Information tab, those sections will be moved to the Baseline Information tab.

Important: the 'Delete' button will delete the patient from the database, including all clinical history related to that patient.

Contact Information

Baseline Information Data Entry - - 039001000000149

Baseline Information | **Contact Information**

Primary Address

Street:

City: State: Zip Code:

County or Area:

Daytime Phone: Nighttime Phone:

Secondary Address

Street:

City: State: Zip Code:

Current Status

Current Status: Date of death: (mm/dd/yyyy)

New Delete Breast Cervical Save Close

The Contact Information section contains address information for the patient.

The Current Status field is used in the Patient Reminder module. Only "Active" patients will receive patient reminders. (Note: if Current Status is left BLANK, patient is considered 'Active' for the reminder system.) Please see Patient Reminders for more information.

When you are finished entering the Baseline information, you can use the buttons on the bottom to enter Breast or Cervical information. You can also click the New button to add another baseline record.

If your Program has created user-defined fields, there will be an additional tab called 'User Defined', which will contain these additional fields.

Cycle/Procedure Data Entry

Adding Cycles and Procedures

Once baseline information has been entered for a patient, Breast and/or Cervical Screening information can be entered into the database. Cycles and corresponding Procedures can be accessed from the Baseline Information data entry screen, or directly from the Find Patient Dialog, which is described in the next section.

To access Cycles and Procedures, select either the Breast or Cervical cycle button. The system will display the existing cycles and procedures. There are several sections in the cycle/procedure tabs.

Cycle Patient Information

The cycle patient information section contains cycle information for the patient. An example of this section is shown below.

Cervical Cycle - Mary T. Smith - 039001000000150

Cervical Cycle Information | Cervical Cycle Cont. | Cervical Cycle Cont. | User Defined

Cycle Location: Old Location: Cyc Num:

Eligibility

Income Eligible: Medicare/Medicaid: Ins. Available:

Suppress Reminders: Suppress MDE:

Clinical History

Prior Pap: Was Pap documented:

Date of prior Pap:

Source of referral:

Add Cycle Delete Cycle < > 1 of 1 Save Cancel

The bottom of the screen contains record selectors, which indicate the current cycle "record" number for this patient. Note: these record selector numbers won't necessarily match the internal cycle number. The internal cycle number can be located on the upper right of this screen. This will be helpful in identifying records from the individual "errors" file from the MDE Edit Program. Please see Running the MDE edit check program for more information.

You can also add or delete cycles using the corresponding buttons on the bottom left. **Note:** deleting a cycle will delete the current cycle and all procedures associated with the cycle.

When adding cycles, the patient information screen will be 'auto-filled' with information from previous cycles. The latest Pap or Mammogram will be used to complete the Prior Pap/Mammogram fields. Please verify this information before moving on to the 'cycle procedures' screen. Please also note that the Prior Pap/Mammograms fields will not be correct if you are entering data out of order. For instance, if you already have data entered from 1998 and 1999, and you want to add data from 1996, the prior Pap/Mammogram information will be from 1999. You will need to overwrite the information that is automatically filled in for these fields.

When you are finished entering this information, you can click the next tab or press "Ctrl" and "Tab" to proceed to the next section.

Cycle Procedures

The cycle procedures and disposition sections contain procedure and status information for the current cycle. An example of this section is shown below.

Cervical Cycle - Mary T. Smith - 039001000000150

Cervical Cycle Information | Cervical Cycle Cont. | Cervical Cycle Cont. | User Defined

Procedures performed

Procedure	Date Requested	Location	Old Location	Date Per

Cycle Disposition

Workup planned:

Final Diagnosis:

Other Dx text:

Status of Final Dx:

Date of Final Dx: (mm/dd/yyyy)

Add Cycle | Delete Cycle | < | > | 1 of 1 | Save | Cancel

Enter all of the procedures for the cycle in the grid labeled "Procedures performed". Each of these procedures will be linked to the current cycle for historical and MDE export purposes. Any user-defined fields will also be included in the row.

Administrators: please see the "Layout files for the procedure data entry screens" section in the 'Administration' chapter for information on how to configure the layout for this procedures portion of this screen.

Note: the Patient Reminder module uses the "Recommended FUP" field to trigger annual and short-term follow-up. For those cases requiring "short-term follow-up", the "Short term follow-up months" field is used to track the number of months (from the current procedure) the patient should return. Please see the Patient Reminders for more information.

The cycle disposition section contains important cycle status information, including the Workup Planned field, which indicates whether further diagnostic procedures are expected. **This field must be completed for**

all cycles, including those with just a Mam, Pap, or CBE.

Note: Pressing the "F2" key will move the user from the "Procedures performed" section to the "Workup planned" field.

When finished entering procedures and cycle disposition information, click Save to post the changes to the database. If you attempt to close without saving changes, the system will prompt you to save the changes.

Cycle Treatment Information

The cycle treatment information section contains treatment information for the current cycle. An example of this section is shown below.

The screenshot shows a software window titled "Cervical Cycle - Mary T. Smith - 039001000000150". It has four tabs: "Cervical Cycle Information", "Cervical Cycle Cont.", "Cervical Cycle Cont.", and "User Defined". The "Cervical Cycle Information" tab is active. Inside this tab, there is a section titled "Stage / Treatment Information". This section contains three input fields: "Stage of tumor:" with a dropdown menu, "Treatment status:" with a dropdown menu, and "Date of treatment:" with a date picker (mm/dd/yyyy). Below these fields is a large empty text area. At the bottom of the window, there is a toolbar with buttons: "Add Cycle", "Delete Cycle", navigation arrows, "1 of 1", more navigation arrows, "Save", and "Cancel".


Please refer to the MDE Data User's Manual for the skip patterns for "Treatment status". For example, this field **must** be completed for certain final diagnosis results, but **may** be completed for other results if treatment was planned.

Note: The User defined tab will contain any fields created by the administrator using the User-defined fields setup function.

Find Patient/Patient Status

Find Patient Dialog

To find patients already in the database, select "Find Patient" from the main menu or from the Search menu bar item. The following dialog will appear.



The "Find Patient" dialog box features a "Field" section with five radio buttons: "Patient Name" (selected), "Social Security Number", "Patient Id", "Chronic Disease #", and "Date of Birth". To the right is a "Value" text input field. Below the field selection is a "More >>" button. A "Search / Refresh" button is located to the right of the "Value" field. The results are shown in a table with columns: "Patient_Id", "Name_Last", "Name_First", "Name_Middle", and "Name_Maia". The table is currently empty. At the bottom, there are buttons for "Status", "New Pat.", "Baseline", "Cerv. Cycle", "Brst. Cycle", and "Close".

You can search by patient name, Social Security Number, Patient ID, Chronic Disease number, or Date of Birth. Select the field on which you wish to search, enter what you are looking for in the Value box, and select the "Search/Refresh" button, or hit the <Enter> key. Results will be displayed, including patient ID and patient name.

If selecting by name, you can search by last name AND first name by separating the names with a comma. For example, to search for Mary Smith, you could enter "Smith,Mary" or "Smith,M" or "Smi,Mary" or even "S,M".

Except for Social Security number and Date of Birth, the system will search for the specified value using "Begins with" logic. In other words, if the field begins with the specified value, the record will be displayed. For example, if you are searching for a Patient Name and enter "Johns" for the value, any name that begins with "Johns" will be displayed.

By default, the results are displayed by either Patient ID or Last Name order, depending on the search. To sort the results by other fields, click

on the field name. You will notice a triangle to the right of the field name, indicating either ascending or descending order. Clicking the field name again will either remove the triangle (no sort), or make the sort order descending or ascending. After making this selection, click the search/refresh button to display the results in the specified sort order.

Once displayed, you can access the data entry screens for a particular patient by selecting the patient (designated by the triangle to the left of the row) and pressing the Baseline, Brst Cycle, or Cerv Cycle buttons.

You can also access the Patient Status screen by selecting the Status button. The patient status screen displays a complete history of the patient, and also provides access to the data entry screens.

If the patient is not in the database, you can also click on the 'New Pat.' button to add a new patient.

Patient Status

Patient Status

The Patient Status screen provides a comprehensive view of a patient's data, using the drill down, or 'branch' method of displaying data. The left-hand side of the screen is used to navigate the different areas, or branches, of a patient's data history. The right hand side is used to display detailed information related to what is selected on the left.

The first 'branch' on the left-hand side is 'Current Status'. This provides a brief overview of the patient, including cycle summary information and last procedure performed. An example is provided below.

Patient Status - 039001000000150 - Mary Smith

<p>Current Status</p> <ul style="list-style-type: none"> Baseline Cervical <ul style="list-style-type: none"> Cycle #1 <ul style="list-style-type: none"> Pap test 1991 - 05/06/2001 - (2) Pelvic Exam - 05/06/2001 - (3) Colposcopy with biopsy - 05/15/2001 Breast <ul style="list-style-type: none"> Cycle #1 <ul style="list-style-type: none"> Clinical Breast Exam (CBE) - 02/03/2000 Mammogram (initial) - 03/04/2000 Cycle #2 <ul style="list-style-type: none"> Additional mam views - 08/09/2000 Biopsy - 08/09/2000 - (2) 	<p>Patient Information</p> <p>Full Name: Smith, Mary T SSN: 000-00-0000 Location: Site 1 Date: Comments:</p> <hr/> <table border="0"> <tr> <td style="text-align: center;"><u>Cervical Summary</u></td> <td style="text-align: center;"><u>Breast Summary</u></td> </tr> <tr> <td>Date of Last Pelvic: 05/06/2001</td> <td>Date of Last CBE: 02/03/2000</td> </tr> <tr> <td>Result: Normal</td> <td>Result: Discrete palp mass - Susp for Ca</td> </tr> <tr> <td>Date of Last Pap: 05/06/2001</td> <td>Date of Last Mammogram: 08/09/2000</td> </tr> <tr> <td>Result: Low grade SIL/HPV</td> <td>Result: Not needed</td> </tr> <tr> <td>Last Procedure: Colposcopy with biopsy</td> <td>Last Procedure: Biopsy</td> </tr> <tr> <td>Date: 05/15/2001</td> <td>Date: 08/09/2000</td> </tr> <tr> <td>Result: Negative (w/NL)</td> <td>Result: Normal breast tissue</td> </tr> </table>	<u>Cervical Summary</u>	<u>Breast Summary</u>	Date of Last Pelvic: 05/06/2001	Date of Last CBE: 02/03/2000	Result: Normal	Result: Discrete palp mass - Susp for Ca	Date of Last Pap: 05/06/2001	Date of Last Mammogram: 08/09/2000	Result: Low grade SIL/HPV	Result: Not needed	Last Procedure: Colposcopy with biopsy	Last Procedure: Biopsy	Date: 05/15/2001	Date: 08/09/2000	Result: Negative (w/NL)	Result: Normal breast tissue
<u>Cervical Summary</u>	<u>Breast Summary</u>																
Date of Last Pelvic: 05/06/2001	Date of Last CBE: 02/03/2000																
Result: Normal	Result: Discrete palp mass - Susp for Ca																
Date of Last Pap: 05/06/2001	Date of Last Mammogram: 08/09/2000																
Result: Low grade SIL/HPV	Result: Not needed																
Last Procedure: Colposcopy with biopsy	Last Procedure: Biopsy																
Date: 05/15/2001	Date: 08/09/2000																
Result: Negative (w/NL)	Result: Normal breast tissue																

View Edit New Delete Print Close


The Patient Status screen also contains summary information for the Baseline data, along with Cervical and Breast cycle data. By selecting each of these 'branches', summary information will be displayed on the right. Information for each 'cycle' and 'procedure' can also be displayed by selecting the corresponding cycle/procedure. For those screens that display procedures, you can also double click on a row to 'jump' to that procedure.

Working with the Patient Status Screen

The Patient Status Screen provides buttons to access the Baseline, Cervical and Breast data entry screens. To edit information for a patient, select the 'branch', and click on the 'Edit' button at the bottom of the screen. For those users with read-only access, the 'View' button can be used to review information.

The 'Print' button will display a summary report for the patient. This report is the same "Patient History with Comments" report found under Standard Reports.

This screen also provides the user with ability to move or copy cycles/procedures. To access these options, select a 'branch' (either a cycle or procedure), and click the right mouse button. The following menu will appear.



The View, Edit, New, and Delete options perform the same function as the buttons at the bottom of the screen. The Cut, Copy and Paste buttons can be used to copy or move information from one cycle to another. For instance, if a procedure was entered into the wrong cycle, you can 'cut' the procedure from the one cycle into the correct cycle. To do this, select the procedure, use the right mouse button, and select 'cut' from the menu. Use the mouse and select the correct cycle. Use the right mouse button again, and select 'Paste'.

You can also perform this same function by dragging and dropping the procedure. Instead of using the right mouse button, select the procedure with the left mouse button, continue to hold the button down and 'drag' the procedure into the correct cycle.

These functions also work between two Patient Status screens. If you identify duplicate patients in your database, open both patient records using the Patient Status screen. Move and resize the screens for easy navigation. You can either use the menu items by right clicking, or use the drag and drop method to consolidate cycles/procedures.

Note: the system will only allow you to copy/paste procedures into the correct 'node'. For instance, you can only copy/move Breast procedures into the Breast cycle 'node'.

Queries

Introduction to Queries

One of the most important aspects of a database system is the ability to retrieve information based on certain criteria. The query or Search Editor


in CaST provides a powerful and easy to use interface for retrieving specific information from any of the tables in CaST.

CaST contains a powerful relational search engine, which supports full Boolean logic. This engine performs the "work" of sifting through all patient records in the database and selecting only the ones that satisfy your specified criteria. In addition, the search editor allows you to select the fields you wish to view and to choose the order in which you want the records to be displayed.

CaST also allows you to save queries for future use. In addition, CaST contains several pre-defined system queries.

Retrieving a Saved Query

To execute a saved query, select Queries from the main menu. The Search Editor window will appear.




Three types of saved queries are available for display: "system" queries, which include standard audits, "all" queries, including those saved by

other users, and queries saved by the current user. You can access these types by selecting the appropriate radio button on the top of the screen.

To execute a query, select the query and select Go. You can also select a query and double click.

Working with Query Results

The following is an example of the results displayed after execution of a query.



The screenshot shows a window titled "Query Result" with a table of data. The table has five columns: Patient Id, Last Name, First Name, Date performed, and Enrollment Location. The first row of data is highlighted in blue. Below the table, it says "1 of 1". At the bottom of the window are buttons for "Export", "View", "Print", "Print Setup", and "Close".

Patient Id	Last Name	First Name	Date performed	Enrollment Location
039001000000151	Valdez	Rosa	11/22/2002	Site 2

The query results will be displayed using the FIELDS and SORT options. The number of records is displayed at the bottom. You can scroll through the results, using the horizontal and/or vertical scroll bars. These scroll bars are only visible if the number of records or number of fields is greater than the viewing area. You can also lock or freeze columns for easier viewing, by right clicking on a field, and selecting 'fix columns'.


You have three other options regarding the results of the query. You can **view** particular patients displayed in the results, **print** the results, or **export** the results.

To access the patient status screen for a patient displayed in the results,


select the patient and select the **View** button. **Note:** This option is only available if the Patient ID field is selected as one of the fields in the query.

To print the results, select the **Print** button. You can also setup several printer options by selecting Print Setup. The following screens will be available.

The image displays two side-by-side 'Print Settings' dialog boxes, both with tabs for 'General', 'Headers / Footers', and 'Margins'. The left dialog box is currently on the 'General' tab and includes sections for 'Borders' (with 'Border' set to '(none)' and 'Border Style' set to 'Solid'), 'General' (with 'Fit to Page' set to 'Never', and checkboxes for 'Repeat Fixed Columns', 'Repeat Fixed Rows', 'Print graphics', 'AutoSize Columns', and 'Center on Page'), and 'Fonts' (with three font selection buttons: 'Select the Table Font', 'Select the Header Font', and 'Select the Footer Font'). The right dialog box is also on the 'General' tab and includes sections for 'Date / Time' (with 'Date Position' set to '(none)', 'Date Format' set to 'dd/mm/yyyy', and 'Time Position' set to '(none)'), 'Title' (with 'Title Position' set to '(none)' and 'Title Text' set to 'CBEs with Discrete Palpable I'), and 'Page Numbers' (with 'Pages Position' set to '(none)', and fields for 'Prefix', 'Separator' (set to '/'), and 'Suffix'). Both dialog boxes have 'OK' and 'Cancel' buttons at the bottom.



To export the results, select the **Export** button. This will open the Save Query dialog, as shown below.



This will allow you to save the results in either a comma-delimited format, or as an Excel spreadsheet. You can use the browse button to select a directory and file name.

Creating a New Query

Creating a New Query -- Getting Started

The first step in creating a new query is determining whether you need to start from scratch. You can often build upon other users queries, even the system queries. If you only need additional fields displayed, a different sort order, or a slightly different subset, you can modify an existing query and save it under a different query name. You can find more information on how to modify an existing query later in this section.

To create an entirely new query, select Queries from the main menu. Select "New - User defined Query" and press "Next".

Selection Criteria

The next step in building a query is defining the selection criteria.


The screenshot shows the 'Query Editor' window with the 'Selection Criteria' tab selected. The window is divided into several sections:

- Fields:** A list of fields including 'B_Workup_Planned', 'Patient_Id', 'BreastCycle_User', 'BreastProcedures', 'BreastProcedures_User', 'CervicalCycle', 'CervicalCycle_User', and 'CervicalProcedures'. 'B_Workup_Planned' is currently selected.
- Operator:** A list of operators including 'is = to', 'is < than', 'is > than', 'contains', 'is <= to', 'is >= to', 'is not = to' (which is highlighted), and 'begins with'.
- Value:** A text box containing 'Planned'. Below it, a list of values is shown: 'Planned', 'Not planned', 'Not yet determined', and 'Planned'.
- Logic:** Buttons for 'And', 'Or', 'Not', and parentheses '(...)' are present.
- SQL View:** A text area showing the generated SQL query:


```
{BreastProcedures.BS_Name} = "Clinical Breast Exam (CBE)"  
AND {BreastProcedures.BS_Result} = "Discrete palp mass - Susp for Ca"  
AND {BreastCycle.B_Workup_Planned} != Planned
```
- Report Title:** A text box containing 'CBEs with Discrete Palpable Mass - susp for cancer - and no workup planned'.
- Owner:** A text box containing 'System'.
- Buttons:** 'Clear All', '<< Back', 'Next >>', 'Go', and 'Exit'.

There are five sections in the Selection Criteria tab: Fields, Operator,


Value, Command Buttons, and Search Statements. The first four sections will determine what is created/displayed in the Search Statements section. The four sections described below will allow you to build one or more search statements:



The Fields section contains all tables and corresponding fields in the database, including user-defined tables. Clicking the '+' to the left of the table name will expand the list of fields contained in the table. You will build your search statement by selecting one field at a time. In the example, the criteria for the first search statement is using 'CBE' procedures. This is a breast procedure, so you would select the BreastProcedures table and select "BS_NAME", which holds the procedure name.




The next part of the search statement is the operator. Using the mouse, select the operator you wish to use with the selected field. In our example, select "is = to" for the operator.



The next part of the search statement is the value you wish to search for

in the selected field. The value you specify is dependent on the type of field. For text fields, like Last Name, you would enter the text. For date fields, you would enter the entire date (in MM/DD/YYYY format – including the '/'). For formatted fields, like our example, the available values will be displayed under the Value box. If you would like to search to see if a field is blank, enter the word **Null** in the value box. When you select the value, click on Add or press the <Enter> key, this will add the entire statement (field – operator – value) to the search statements.

Command Buttons



If you are only creating one search statement, you do not need to use the command buttons. Command buttons are used when you use two or more statements. You can select the buttons before or after you create the statement.


The first three command buttons (AND, OR, NOT) control the 'logic' of the statements. Placing the 'AND' between two query statements will result in records where BOTH statements are true. Placing the 'OR' between two query statements will result in records where EITHER statement is true. The 'NOT' can be used to select the opposite of a statement. For instance, if you are searching for all Races BESIDES 'White', you can build a statement NOT 'Race = White'.

The '(' and ')' command buttons are used to group logic statements together. For instance, if you wanted to search for race equal to White OR Black, along with Ethnicity equal to Hispanic Origin, you would need to use these parentheses to group the two search statements on race. So, the statement would be: (Race1 = White or Race1 = Black) AND (Ethnicity = Hispanic Origin). You can add or remove these parentheses at any time. To remove, right click on the search line, and select "remove".

The "Clear All" button at the bottom clears all search statements. Selecting a statement and using the right mouse button also accesses the delete and insert line functions.

Helpful hint: for complicated queries, it is sometimes easier to first create all your search lines, then go back and add the 'commands' (AND, OR, NOT), along with the right and left parentheses.

Fields Selection



The Fields Selection tab allows you to specify what fields are displayed in your query. All fields in the database are listed on the left and the selected fields are listed on the right. The buttons in the center are used to select and de-select fields. The arrows indicate which direction the field(s) will move. The double arrows are used to move ALL fields from one side to the other.


The "Move Up" and "Move Down" buttons are used to order the selected fields. To move a field, select it in the right hand column and click on one of the move buttons.

IMPORTANT: in order to access the 'Patient Status' screen from the Query Results window, 'Baseline.Patient_Id' must be selected as one of your fields. CaST will automatically select this field when you create a new query.

Note: the number of records displayed is partially determined by which fields you choose on this screen. For example, if you are searching for a baseline field (i.e. Race), but choose fields from the Breast Procedures


tables, you will most likely get multiple records per women, since the procedure tables hold all procedure records for the woman.

Sort Options



The Sort tab allows you to control the sort order for the displayed records. The "Order by" box must be checked for the selected sort order to take affect. The fields selected in the previous Fields tab will be available on the left-hand side of the tab. Moving fields from one side to the another is accomplished in the same way as in the previous tab.

SQL Tab



This final tab allows you to view the SQL code which constitutes the query and save the query for later use. You can also specify a report title and select Go from this tab. Note: you can select a title and Run the query from any of the tabs in the Search Editor. When you select Go, the results of the query will be available for viewing, printing and saving. Please see Working with query results.

The Edit SQL button is used to modify the SQL code. This is explained in the Advanced SQL section.

Modifying an Existing Query

If you have an existing query, either a user-defined or system query, which you would like to modify, CaST allows you to modify the query and save it as one of your user-defined queries.

To modify a query, select the query from the Search Editor screen. Scroll through the various tabs and make changes as necessary. In the SQL tab, click save. If you are the owner of the query (i.e. you originally

created the query), you can overwrite the existing version. If you are not the owner, you must specify a different name for the query. Queries that you create/save will show up when you select the "User.." button on the Search Editor screen.

Note: the 'owner' of a query will also appear at the bottom of the 'Queries' tab, as you select a query.

You will notice that some queries do not have the normal search statements, field, or sort options available for editing. These queries were created using Advanced SQL. You can still modify these queries, but you need to do all of the modifications in the SQL tab, which requires a familiarity with SQL code.

Advanced SQL

The SQL tab in the Search Editor contains an Edit SQL button on the top left hand side of the screen. This button can be used to edit the SQL code for any query. To use this feature, you must have an understanding of SQL or have access to an SQL editor, such as Microsoft Access 97. If you use Access 97 to build your own queries on the data tables, you can cut and paste SQL code into the SQL tab. This is necessary to execute complex queries such as the System Query "Calc: Average Number of Days: Breast Procedures". This query uses the Group By option and calculates Min and Max values.

When using this feature, please be aware that the Selection, Fields and Sort tabs will not be available for future modifications. Please also be aware that selection criteria must be created using valid SQL code.

Note: when creating criteria statements, you must use single quotes around values instead of double quotes (by default, Access uses double quotes for 'text' fields).

Standard Reports

What are Standard Reports?

Standard Reports are an extension of the Queries function in CaST. The difference between queries and standard reports is that standard reports use a graphical reporting package and are displayed in an easier to read format.

The standard reports window is very similar to the queries or Search Editor window. The differences are that you can't create a new standard report from scratch and you are limited to predefined tables and fields in the selection criteria tab.

There are several "system" Reports provided with the installation of CaST. You can run these reports directly or you can run the reports using a different subset, by modifying the Selection Criteria.

Running a Standard Report

To run a standard report, select 'Standard Report' from the main menu. The Report Editor window will appear.

Three types of saved reports are available for display: system reports, which include those provided with the installation, all reports, including those saved by other users, and reports saved by the current user. You can access these types by selecting the appropriate button on the top of the screen.

To execute a report, select the report and select Go (or double click). A basic report window is displayed. You can view the results and print the output from this window.

Another option available when running a report is the ability to create a PDF file. To create a PDF file instead of viewing the report in the standard window, check the 'create PDF file' box and select Go. You will be asked to provide the name and location for the new file. If you have Adobe Acrobat on your machine, you can then open the PDF file.

Modifying an Existing Standard Report

Although you can't modify the design of a standard report, you can modify the underlying subset or sort order. You can use the Search Editor screen the same way you would to modify Queries. You will not be able to change the tables or fields related to the report.

After making the changes to the Selection Criteria or Sort tabs, you can save the options by renaming the report in the SQL tab. The current report with your modified changes will be saved for future use.


Patient Reminders

What are Patient Reminders?

The Patient Reminders module allows the user to generate a list and corresponding mailing labels for all women due in for services. Users can generate reminders for Breast and/or Cervical screening services.

To successfully run the Patient Reminders module, the "Recommended Follow-up" field must be completed.

How to Run Patient Reminders

The screenshot shows a Windows-style dialog box titled "Patient Reminders". It is divided into two main sections. The left section, titled "Screening", contains two checked checkboxes for "Breast" and "Cervical". Below these is a "Time Interval" section with "Starting Date:" and "Ending Date:" fields, each followed by a date input field (MM/DD/YY). Underneath is an "Enrollment Site" section with a dropdown menu currently set to "All Sites". The right section, titled "Type of Reminder", contains three checked checkboxes: "Annual", "Short-term Followup", and "First Screening". Below "Short-term Followup" is a numeric input field with "6" and the text "Default # of Months". Below "First Screening" are two more numeric input fields: "12" with "Months from Enrollment" and "50" with "Mammogram (age)". Below these is another numeric input field with "40" and "Pap Test (age)". At the bottom right of the dialog is a checkbox for "Create PDF Report", which is currently unchecked. At the very bottom are three buttons: "View", "Run" (which has a mouse cursor hovering over it), and "Close".

Patient Reminders can be run for Breast and/or Cervical screenings. Click the appropriate **Screening** type boxes to generate lists and labels.

The starting and ending dates for the **Time Interval** section indicates the range for when the women are due for services. For example, to generate reminders for women due in for services during the month of March 2003, you would enter "03/01/2003" and "03/31/2003" in the Starting and Ending Date fields, respectively.

Choose "All Sites" or the individual **Enrollment Site** in which you want to generate reminders. Note: a future version of the system will allow you to choose multiple enrollment sites.

You can generate one or more of the following **Types of Reminders**: annual, short-term follow-up, and first time screenings.

Annual: patients with "follow routine screening" indicated

Short-term: patients with "short-term follow-up" indicated. If the "short-term FU" months field is not completed in the database, the default number of months specified on this screen will be used in the calculation.

First screening: patients with no screening visits performed

The "Months from enrollment" field is used to ensure that a "first screening" reminder is not sent to women who may have already had a screening, but the information is not in the system. This number can be adjusted according to your programs data entry schedule and how enrollments are performed.

The Mammogram and Pap Test Age boxes are used to indicate the minimum age to warrant a first time reminder. For example, if Mammogram Age = 50, women who turn 50 during the interval will have a reminder generated. Also, for other first time screenings, only those women over 50 will receive a reminder.

Click on "Run" to generate reminders. A message will be displayed when the reminders are processed. Click "Ok". Two windows will then automatically open. One is a report containing the mailing labels, the other is a standard "Query Results" window containing information on the women.

Note: as in standard reports, you have the option to send the output to a PDF file. To do this, check the 'Create PDF File' box.

Working with Patient Reminder Results

The Patient Reminders module creates a temporary table on the user's hard drive. This table is cleared each time the module is run. After the process is complete, a query results window and the mailing labels will appear.

The mailing labels will appear as a report and can be printed using Avery 5160 stock labels. Each label will contain coding to represent the reason for the reminder letter. Each label can have a breast and/or cervical reminder code. The coding scheme is as follows:

ba: breast annual
 bs: breast short-term follow-up
 bff: breast first screening (over the minimum age, i.e. ≥ 50)
 bfa: breast first screening (birthday during the interval period, i.e. 50th)

 ca: cervical annual
 cs: cervical short-term follow-up
 cff: cervical first screening (over the minimum age)
 cfa: cervical first screening (birthday during the interval period)

The query results window will contain the same women identified on the mailing labels. Identifying information, including Patient ID, name and address will be included in this table/query results window. The other fields of interest include the baseline date, performed date, breast reminder and cervical reminder fields. These fields are described below.

Baseline date: left blank for first screenings. For short-term and annual reminders, screening exam date is used. If both types of screening reminders are due, two baseline dates (separated by a ‘;’) are specified.

Performed date: only completed for short-term follow-up.
 Performed date for procedure that called for short-term follow-up is used to populate this field.

Breast/Cervical Reminder Fields

The BreastReminder and CervicalReminder formats in the cast.fmt library are used to populate these fields. The formats are as follows:

1=1st screening-age (patient has birthday during interval period)
 2=1st screening (patient is due for first screening)
 3=Short-term follow-up:
 4=Annual follow-up

The results can be viewed, printed, and exported similar to any other query. Please see Working with Query Results for more information.

Note: previous versions of CaST stored the reminders results in a table in the main database. CaST now uses temporary databases for this function. You can, though, export the reminder results into either an Access97 or Access 2000 database.

How Patient Reminders are Created

The patient reminders module uses several pieces of information to determine if a patient is marked for inclusion. The following describes the fields and logic in determining which patients are included.

Current Status Field (contact information tab on Baseline record)

This field must be either set to "Active" or left blank.

Suppress Reminders (cycle information screen)

If this field is set to "Yes", patient reminders will not be generated for that patient AND cycle. It is possible to have one cycle set to "Yes" and the next set to "No".

First screenings

Women with no procedures performed who meet the minimum age requirements will be included. The women must be in the system (using the enrollment date and interval start and end dates) for at least the number of months indicated on the Patient Reminders screen. Note: date of birth and enrollment date must both be completed for these types of reminders to be generated properly.

For example, if the minimum age is "50" for Mammograms, and the number of months is set to "12", the following cases will be included.

- Those women who have been enrolled for at least 12 months and have their 50th birthday in the interval.
- Those women who have been enrolled for exactly 12 months and are at least 50 years old.

Annual screenings and short-term follow-up

The reminder module searches for the latest cycle, using the performed dates. It then searches for the latest procedure performed in that cycle. It will then look at the "Recommended FUP" field for that procedure. If two procedures are performed on the same day, a short-term follow-up will over-ride a "follow routine screening".

For "Follow routine screening" recommendations, the performed date for initial screening exam is used to calculate the 12 month rescreen date.

For short-term follow up, the procedure date and the number of

months specified in the "short-term FU months" fields are used to calculate the date for the follow-up. If the number of months is not specified in the database, the default number of months specified on the Patient reminders screen is used.

Other Reports

Overdue Reports

What is an Overdue Report?

The Overdue Report is a feature taken from previous versions of CaST. This report allows the user to query the database for "overdue" information, including results, notification, and follow-up data.

The Overdue Report dialog allows the user to select the type of overdue information and number of days for both breast and cervical screening information.

Note: There are also system queries, which pull similar information from the database. You will have to run separate reports and will not be able to specify the number of days, but you will be able to quickly access the patient records from the query results window.

Please also note that some of items in the overdue report use optional fields. For instance, the 'notify client' section uses the 'notify date' field, which is NOT a required field. If your Program does complete this field, a report with this option will include all records in your database. Also, the 'date requested' field is queried in the 'results' option for overdue items. 'Date requested' is also not a required field.

How to Run the Overdue Report

Overdue Report

Type of Report

- ☒ Breast Screening
- ☒ Cervical Screening

Select a Site

All Sites

Overdue Items

- ☒ Results 60 Days
- Notify Client**
 - ☒ Normal Screening 30 Days
 - ☒ All Others 30 Days
- Cycle Information**
 - ☒ Diagnosis 30 Days
 - ☒ Treatment 30 Days

Group By

- ☐ Site

☐ Create PDF Report

☐ Override Report Title

Report Title: Breast Screening and Cervical Screening Overdue Report - All Sites

Run Close

To run the overdue report, you must select the options you wish to include in the report. You can select Breast and/or Cervical information, along with the options of Results, Notify Client and Cycle Information. You can also specify the number of days used in determining if an item is "Overdue". You can also run the report for all enrollment sites and sort the results by site. The report title is automatically created, but you can override the title by checking the "Override Report Title" box and supplying the title.

Note: as in standard reports, you have the option to send the output to a PDF file. To do this, check the 'Create PDF File' box.


Duplicate Patient Report

Finding Duplicates in the Database

CaST contains a utility to search for duplicate patients in the database. This utility uses logic that looks at several variables, including Social Security Number, name, and birth date to determine if duplicates exist.

Running the Duplicate Patient Report

To run the duplicate patient report, select SEARCH|Find Duplicate Patients from the Main Menu. The following screen will appear.




The 'Min percent match' field is used to adjust the sensitivity of the 'Find Duplicate Patient' function. This function loops through the database, looking for matches in several fields between patients. Setting this 'percent match' to 60 indicates that you want the function to return all records which have a 60% match, i.e. 3 out of 5 key fields the same. The higher the number, the less possible duplicates will be reported.

To start the check, click on Run. The function will make four passes through the database. When the process is complete, you can view any duplicates found by clicking on "View Results". This will bring up a screen similar to the query results window, but will group any duplicates using color coding. If duplicates are found, you can then click on View and use the Patient Status screens to consolidate records.

MDE Export

Creating the MDE Export File

To create the MDE export file, click on the Generate MDE button on the main menu, or choose Tools|Generate MDE from the menu bar. The following screen will appear.



The "Submission Date" and "Include data through" parameters will be automatically completed according to the current date. The system will calculate the most logical submission date and adjust the "Include data through" accordingly. Please confirm these dates using the information in the "Key Dates" section in the Data User's Manual.

The "Use Patient Names..." will create an MDE export file with patient name information included in the ID field, instead of the system Patient_Id. This feature can be helpful when running the MDE edit program, but should **NOT** be used when creating the final file to be submitted to CDC.


To create the MDE file, click on Export. A status bar will indicate the current progress of the process. The number of records written to the MDE file will be displayed with the process is complete.

When completed, you will be able to view the file by clicking View File. This will bring up a simple text viewer with the MDE file.

To run the MDE edit check program, click on the "Run Edit Check" tab.

Running the MDE Edit Check Program

After you Create the MDE export file, the three parameters for the MDE file, Summary File, and Error File will be initialized in the Run Edit Check tab.



To initialize other MDE edit check parameters, you should click on "Edit INI file" and bring up the INI file to edit. Please refer to the MDE Edit Check documentation for information on the edit check program and the INI file.

NOTE: You may need to edit the Program_Abbreviation entry in the INI File before the first time you run the Edit Check program.

To run the edit report on one edit message/number, check the box and specify the edit type and number. To bypass prompts in the edit program, check the "Bypass prompts" box.

To run the edit check program, click on the Run Check button. This will launch the edit program in a separate window. When the edit check program is completed, you can close the window.

To view the summary or individual edit/error file, click on the appropriate button. The MDE file, along with the summary and error log files can be found in the MDE folder off of the Cast II directory.

The 'View Errors' button will allow you to view the 'individual edit file', which contains the patient ID information. With the 'Generate MDE' window open, you will not be able to access other areas of the system, including the 'Find Patient' screen. You can, though, leave the 'individual edit file' window open and close the 'Generate MDE' window. This will allow you to open the 'Find Patient' window and 'copy/paste' patient IDs between the two windows.

Note: when viewing the individual error listing, please note that the record ID field is created using the following scheme:


1st digit: 1=Breast, 2=Cervical

Last 5 digits: created using the cycle number field in the corresponding cycle table. This internal cycle number may not match the record locator number on the bottom of the cycle data entry screens. The cycle number field can be located on the top right of the cycle data entry screen.


IMPORTANT: after you use the reports of the MDE edit check program to correct any errors, you must first re-create the MDE export file before you re-run the edit check program. The edit check program uses the MDE export file, so if you don't first re-run the export, you will not see any difference in the error reports.

Compressing the MDE File for Submission

The MDE text file should be compressed in preparation for submission. To compress the MDE file, click on the "Compress MDE File" tab. The following screen will appear:



Click on the "Compress". This compression utility will zip the MDE file (.txd) and create a ".gz" file. When completed, a message similar to the following will appear:



Please note the location and name of the compressed file. You will need to either copy this file to a diskette or attach it to an email for submission.

Working with Decentralized Data Entry

Import/Export Functions

What are the Import/Export Functions?

For programs utilizing decentralized data entry, CaST uses the Export and Import functions to transfer data between the decentralized data entry sites and the central site. These functions can be accessed from the Tools menu.

Export Data to the Central Site

To Export (or Send) data to the central site, select Export from the Tools menu bar off of the main menu. This will bring up the following screen.

The screenshot shows the 'Export' dialog box with the following fields and options:

- Export Subset Selection:**
 - ☒ Enrollment Sites: All (dropdown menu)
 - ☐ Patient Id: (empty text field)
 - ☐ Date Range: Start Date (empty date field) 00:00:00, End Date (10/24/2002) 23:59:59
- Output Database:**
 - Database Type: Access97 (dropdown menu)
 - Name: C:\CaST Files\Export\001_All_ES_20021024.mdb (text field with browse button)
 - ☐ Zip Export File
- Buttons:** Print Log, Export, Close

A decentralized site can export by enrollment site, patient ID or date range. Choosing "All" for enrollment site will export all data in the current

database.

The default export database name will be created as you choose the export specifications, and will be created using the following algorithm (note: all files will begin with the corresponding Data Entry Site code – i.e. 101, 102, etc):

Enrollment Site: Enrollment Site (or "All") + "_ES_" + current date
Patient ID: Patient ID + "_PI_" + current date
Date Range: Start Date + "_" + End Date

You can keep the database name and location or modify it as needed. Click on "Zip Export File" if you wish to compress the file while you are exporting the data. This is helpful if you need to send the database using diskette or email.

Note: you should **not** export directly to a diskette, as the file exported may not fit on a diskette. You should export to a larger media device, and then copy the export file or zip file to the diskette, if it is small enough to fit.

Click on Export to begin the export. When completed, the 'Export Database' file can then be sent to the central site for inclusion in the master database.


Note: the messages that appear on the second 'tab' of the screen are also written to the 'export log', which can be accessed at the specified location or printed using the 'Print Log' button. These messages include the number of records exported from each table.

IMPORTANT: you will not be able to export records that do not 'belong' to your decentralized site. If a woman transfers from another site and you used the import function to transfer one or more patients, you will first need to use the 'convert patient IDs' function, before you export your data.

Import Data from a Data Entry Site

Each decentralized site will send the central site an export database. Since these files MAY be named the same, it is up to the database administrator to rename or organize these files to avoid data loss. Note: if zipped files are received on diskette, they should first be moved to the hard drive or network.

First, make sure no one is using the system. To run the import, select Import from the Tools menu. The following screen will appear.



You can select "MS Access Database (*.mdb)" or "Zip Files (*.gz)" to search for files. Select the database or zip file(s) to process and click on Import. Note: to select more than one file, hold the 'CTRL' key down and select the files.

The IMPORT function will first execute several checks on the file(s) selected to be imported. Before any data is imported, the user will have the opportunity to view any errors, using the following screen:

Import Verifications

Please review each Import Verification. Anything in red would prevent you from importing that database. Please view log for more detail.

☒ **C:\Cast Files\Export\101_All_ES_20020916.mdb**

* Table Count Verification for Site 101:

Table Name	Main MDB	Import MDB
Baseline	146	147
Baseline_User	146	147
BreastCycle	147	147
BreastCycle_User	147	147
BreastProcedures	242	243
BreastProcedures_User	242	243
CervicalCycle	146	147
CervicalCycle_User	146	147
CervicalProcedures	141	141
CervicalProcedures_User	141	141

☒ **C:\Cast Files\Export\301_All_ES_20020916.mdb**

* Table Count Verification for Site 301:

Table Name	Main MDB	Import MDB
Baseline	146	146
Baseline_User	146	146
BreastCycle	147	147
BreastCycle_User	147	147
BreastProcedures	242	242
BreastProcedures_User	242	242

Any errors will appear in red, and will prevent that database from being imported into the main database. You can review the information provided on this screen, and determine which files you would like to be imported. This verification screen includes the number of records currently in the database for the site, along with the number of records in the database to be imported.

Note: the messages that appear on the second 'tab' of the screen are also written to the 'import log', which can be accessed at the specified location or printed using the 'Print Log' button.

The provider sites from the import file will be updated in the central database, except for the program wide list of providers, indicated by the "001" prefix. Please see Provider/Enrollment List Setup in the Administration section for more information on program wide providers.

Restoring/Creating New Site Databases

Creating a Restored or New Site Database

The central site can easily restore a site's database, in the event that the site loses their database, or their database becomes corrupt beyond repair. The central site can also create a 'shell', or empty, database for any new sites that begin screening and/or performing data entry.


Note: along with the 'patient' database that would be created in the following functions, the central site would also need to provide the 'users' database (users.mdb) and the 'users.fmt' file (if any user formats were created).

These two functions can be accessed through the Administration menu, by selecting 'Decentralized Setup'.

Restoring a Site's Database

Since the central site will receive updates from each of the decentralized sites, a 'restored' database can be created from the central database. This would be used if the site loses their database and can not restore from a backup.

The first tab on the Decentralized Setup screen is the 'Restore Site DB' function:



Select the site to restore. The available sites are those that were entered into the 'Site Information' window on the 'Administration' menu. Click on 'Restore' to restore the database.


The name of the restored file is listed under 'Database Name' and is comprised of 'castmain_restore', the site number (i.e. 101) and the current date (in YYYYMMDD format). This function will also create a compressed version of the database, with the extension 'gz'.

You will send the 'gz' file to the site. The site would copy the file to the hard drive and then use the 'Restore' function off of the 'Tools' menu to restore their data.

Creating a New Database for a Site

If a new data entry site is created, the central site will need to make sure that this new site is using the same database structure as the rest of the sites. The new data entry site would install CaST, and then use a file created by the central site.

The second tab on the Decentralized Setup screen is the 'Create New/Empty Site DB' function:



The screenshot shows a window titled 'Decentralized Setup' with two tabs: 'Restore Site DB' and 'Create New/Empty Site DB'. The 'Create New/Empty Site DB' tab is active. It contains the following elements:

- Enter Site Number:** A text box containing the number '105'.
- Progress Bar:** A box labeled '0 % Done'.
- Database Name:** A text box showing the path 'C:\Cast Files\Export\castmain_create105.20030116.mdb'.
- Database Type:** A dropdown menu currently set to 'Access97'.
- Zip Empty Database File:** An unchecked checkbox.
- Buttons:** 'Create' and 'Close' buttons at the bottom right.

Enter the site number for the new site. This is the same number that you will assign in the 'Site Information' window on the 'Administration' menu. Click 'Create'.

The name of the database is listed under 'Database Name' and is comprised of 'castmain_create', the site number (i.e. 103) and the current date (in YYYYMMDD format). This function will also create a compressed version of the database, with the extension 'gz'.

You will send the 'gz' file to the site. The site would copy the file to the hard drive and then use the 'Restore' function off of the 'Tools' menu to overwrite the default 'empty' database created with the CaST install.

Exporting/Importing Provider Sites

Using the Provider List Export and Import Functions


The provider list export and import functions will primarily be used to update program/state wide provider lists at programs utilizing decentralized data entry. The central site can maintain this list, and if necessary, export an updated list to each of the local data entry sites. The local site can then import the sites and update the state/program wide provider list.

The export function can also be used by any program to access provider information, such as name and address.

Both functions can be accessed from the Provider List option on the File menu.

Exporting Provider Sites

To export one or more provider sites, select Export from the Provider List submenu off of the File Menu. The following screen will appear.



You can select providers from all sites or those from specific data entry sites. To view/select providers from a specific data entry site, select the data entry site. The available providers will show up on the left-hand side of the screen. To select providers to export, highlight the provider and select the > button, which will move the selected item to the right hand side. You can move multiple providers over by using the >> button. Use the < and << buttons to un-select providers.


To select an export file name, type in the full path name or click on the "..." button. The default file extension will be "TXT". Click on Export.

Two files will be created, using the filename specified. The "TXT" file contains the providers, and a "ZDS" (also a text file), will contain the file layout, or description.

The "TXT" and "ZDS" files can be sent to the local or central site for import into the database.

Importing Provider Sites

To import one or more provider sites, select Import from the Provider List submenu off of the File Menu. The following screen will appear.



Enter the full path and name of the import file (and ZDS file) or click on the "... " button to select the file. If you want to update any providers for your data entry site, select "Update information for current site". If you don't have this checked, provider sites from other sites (i.e. program wide providers, other data entry sites) will be updated.


Click 'Import' to begin.

Administration

User Setup

To set up user accounts, you must be logged in with administrator privileges. The default login (**Admin**) has these privileges. Select "Administration | User Information" from the Main Menu. The following screen will appear.

Note: Programs using the decentralized data entry method should first setup the Data Entry Sites before creating users.



You can edit or delete existing users. To add a new user, click 'New'. The following screen will appear.

Login:
 Name:
 Phone Number:
 User Number: Site Number:
 Site Name:
 Password:
 Verify Password:

Permissions:
☒ Read
☐ Write
☐ Delete
☒ Edit Providers
☐ Administrator account

The following are the fields:

- Login:** Login used to gain access to CaST.
- Name:** Full name of person using the Login.
- Phone Number:** <Optional> Phone number of user.
- User Number:** Unique number to identify the user – system creates.
- Site Number:** Used in Decentralized installations. For centralized installations, use the default 001.
- Password:** Initial password of user (user can change later – either leave blank or make it the same as the user name)
- Verify Password:** Verification of the password entered.
- Permissions:** The permission level of the user.

Edit Providers: Can this user add/modify/delete enrollment sites, primary providers, or providers?

Administrator: Is this user an administrator? If so, he/she will have the ability to create formats, user-defined fields, add/remove users, and have access to the Administrator menu

Save the user by clicking 'Save'. Add additional users in the same manner, as necessary.

Provider/Enrollment List Setup

The ability to edit provider lists is a new feature in CaST. Each data entry site can have a list of providers. These providers fall into three categories: Enrollment Site, Primary Provider, or Provider. Check boxes on the provider list setup screen will determine whether a provider is defined as one or more of these provider types. These provider types will also determine the values that will displayed for the enrollment site, cycle location and procedure location fields, respectively.

To setup the provider, select File | Provider List from the main menu. The following screen will appear.

Site	Name	Status	Contact
0010001	Site 1	Enroll (I), Primary (A), Provider (A)	
0010002	Site 2	Enroll (A), Primary (A), Provider (A)	

By default, all providers are displayed. You can use the filters at the top of the screen to subset the list by type and/or status.

The administrator, or a user with the 'Edit Providers' box checked in their user properties, can edit or delete providers from this screen. To add a new provider, press the New button. The following screen will allow you to enter information.

Enrollment/Provider Properties

Number: 0010001

Name: Site 1

FIN, Site #:

Contact:

Address:

City, State, Zip: Columbus OH

County/Area: Adams

Phone, Fax:

☐ Enrollment Site (Inactive)

☒ Primary Provider

☒ Provider

New Save Close

The provider number will automatically be created using the data entry site and a sequential number. For programs utilizing centralized data entry, all provider numbers will begin with "001". For programs utilizing decentralized data entry, "001" will hold any program/state wide providers. These will be created before the database is rolled out to the data entry sites.

The three check boxes in the upper right corner will determine the type of provider. The "Name" will be used for data entry. When completed, either Save the record, enter a New provider, or Cancel your changes.

It is very important that you complete the City, State and County for each entry. These fields will be used in the creation of the MDE export file.

FIN, or Federal Identification Number, and Site # were added at programs request. These fields are optional.

When changing the 'type' for existing providers, the system will first check to see if records already exist. If records exist for a provider, the provider will be marked as 'inactive'. The provider will still appear in the drop down lists, but will have the tag of '(inactive)' after the name. Also, all 'inactive' providers will appear at the end of the drop down list, so that data entry staff can easily view and use only active providers. The type ahead feature available in CaST will therefore not select these inactive providers.


Note: There is a Re-Sync function on the File | Provider List menu. Most users do not need to worry about this function. This function was added to CaST 5.0 for those programs that have external provider databases. Since CaST 5.0 has additional fields in the provider table to keep track of 'inactive' providers, this function was added to help those programs that

have not yet updated their external provider databases for these new fields. When run, this function automatically checks the entire database and makes sure that these new provider status (inactive versus active) fields are set properly.

Data Entry Sites

If your program is using the decentralized data entry method, you need to setup the data entry sites. Programs using the centralized data entry method do not need to modify this data. The system contains a default data entry site.

Each decentralized site will be provided with a site number, which will be contained in other system tables. The decentralized site will also use this number during installation of the system. As the centralized site, you will assign the numbers given to these decentralized sites, along with address information. To setup these sites you must be logged in with administrator privileges. From the main menu, select **Administration | Site Information**. The following dialog will appear.



You can edit or delete sites from this screen.

Important: the central site should always be site "001". This will allow the central site to create program wide provider sites.

To add a site, click 'New'. You will now be able to enter site information using the following screen.

The screenshot shows a 'Site Properties' window with the following fields and values:

- Site Number: 002
- Site Name: Data Entry Site #2
- Address Line 1: 123 Main Street
- Address Line 2: (empty)
- City, State, Zip: Baltimore MD 99999
- Phone: (555) 555-1212

Buttons at the bottom: New, Save, Cancel.

The following are the fields:

Site Number:	Unique 3 digit code - example 002
Site Name:	Name for the site
Address Line 1:	Address
Address Line 2:	Additional address information
City:	City
State:	State abbreviation
Zip:	Zip code of site
Phone:	Phone number of the site

Save the record by clicking 'Save'. Add additional sites in the same manner, as necessary.

Layout Files for the Procedure Data Entry Screens


The administrator can configure the layout of the procedures data entry tab, or procedures 'grid'. The administrator can control which fields are displayed, the field width, and the order in which the fields are displayed. This layout will be available for all users and can be configured for both the breast and cervical procedure grids.

To modify the layout for the grid, go to the cycle procedures section.

You can change the order in which the fields are displayed by selecting and moving the field column. To do this, point the mouse cursor over the field name until the 'down arrow' appears, and then click the left mouse button. You can then 'drag' the field to the desired location by holding down the left mouse button and moving the field to the correct location.

You can also change the width of the field. To do this, position the mouse over the vertical line that separates the field (to the right of the field name). The cursor will change and you will see one vertical line, with two horizontal arrows. Use the mouse and left button to re-size the field.

You can also hide columns. To access this option, along with the other options discussed below, you will need to right click in the procedure grid. The following menu will appear:



To hide a column, select the column (by pointing the mouse cursor over the field name), right click, and then select 'Hide Column'. If you decide later to display the column, you can use the 'Show Column' option.

Once you are completed with the layout changes, right click and choose 'Save Layout'. Layout files should be saved in the directory determined at system installation and have a 'lyt' extension. By default, the system uses layout files called 'breast.lyt' and 'cervical.lyt'. It is recommended that you use these file names.

It is possible to setup additional layout files, and have some users access different layouts. To do this, setup the layout as noted above, but save the layout files with different names. You will then need to logon to the user's machines as the administrator and choose 'Layout file location Editor' from the 'Administration' menu. Choose the new layout file(s). That user will now access the new layout file. You can also temporarily have the user access a different layout file, but right clicking, and selecting 'Load layout'.

By default, the 'Procedure Name' field is locked. This means that as the user scrolls to the left, the procedure name is always in view. If you wish to display additional fields, you can 'lock' (and then 'unlock') columns. Locking a column will effectively lock that column, and all columns to the left. Locking and unlocking can not be saved to layout files, and 'expire' once the screen is closed.

Database Password/Encrypt Editor

CaST 5.0 includes a function to allow the administrator to password-protect and/or encrypt the patient database.

By default, the patient database (castmain.mdb) is not password-protected. Even though users must have a username and password to gain access to the CaST main menu and data entry screens, a user with access to the directory where the databases are located could open the raw data. For this reason, it is a good idea to password-protect the database.

For additional security, it is possible to encrypt the database. Even if the database were password-protected, it could be opened using a plain text editing program. It would be more difficult, but still possible, to view the information in an unencrypted database. An experienced computer user could discover the password as well.

From the Main Menu, select "Administration | Database Password Editor". The following screen will appear.

The image shows a Windows-style dialog box titled "Database Password Editor". It contains three text input fields labeled "Old Password:", "New Password:", and "Verify Password:". Below these fields is a checkbox labeled "Encrypt the Database". At the bottom right of the dialog are two buttons: "Apply" and "Close".

To add or modify the database password, type the current password in the 'Old Password' field, then type the new password in both the 'New Password' and 'Verify Password' fields. If no password has been set previously, leave the 'Old Password' field blank.

To encrypt the database, check the 'Encrypt the database' box. Note that you **must** specify a password if you are encrypting the database.

Adding a password does not affect the performance of the database. In theory, though, encrypting the database could affect performance. This was an issue in the past; however, since modern networks and workstations are faster and more powerful than before, it is unlikely that the effects will be significant. Administrators should monitor the

performance of the database access after an encryption to be certain.

If you no longer want the database to be encrypted, uncheck the box to decrypt it. You will have to enter at least the old password to do this.

Click 'Apply' and your changes will take effect. Click 'Close' to close the Database Password Editor window.

Note: for the same reasons given above, it is also recommended that you encrypt the 'users' database (users.mdb) which stores user passwords, queries and reports. This database is already password-protected. You can encrypt it by running the 'CaSTEncryptUsers.exe' program located in the 'utilities' directory on the CaST 5.0 CD-ROM.

Format Library

What is a Format Library?

CaST uses format libraries to store Label information for fields in the database. Format libraries make it possible to update the display or label definitions for fields in the database without making structure changes to the data or making changes to the system itself. This also gives the user the flexibility to create formats (drop down lists) for any user-defined fields created in the system.

There are three format libraries used the system, the default format file (cast.fmt), the county formats (county.fmt) and the user-defined formats (user.fmt). The only format library which can be modified by a system administrator is the user-defined library.

The User-defined Format Library

The user-defined library will contain all formats created by the administrator. At the request of programs, the formats for race, referral source, and cpt codes are also contained in this library.


When modifying this library, the administrator should be aware of the following issues.

- If adding values, make sure the value doesn't already exist.
- The length of the field must be large enough to handle the value of the format.

- Do not delete values for the Race, MDE_Translate_Race50 or Referral Source formats.
- Deleting values in the format does NOT change the information in the database.
- There are two formats for CPT codes, one for Breast and one for Cervical. There is some overlap, particularly with the office visits.
- When adding values for the Race format, you **must** also add a corresponding entry in the MDE_Translate_Race50 format. This translation is used for the MDE export program, indicating which MDE category the new format will be placed in. For example, you could add a value of '10' to the Race format, and label it "sample Race". If this race should be reported in the "Asian (3)" MDE category, you would also create a row in the MDE_Translate_Race50 format - '10' in the value column, '3' in the label column. This will tell the MDE export program to report all '10' values to '3' - Asian.
- Using a format in a user-defined field will create a drop down box and will restrict the user from entering any data values other than those specified in the format.

Using the Format Editor

To view/edit format libraries, select Format Library from the Tools menu. The following screen will appear.



The default format dictionary is displayed. This dictionary contains the cast.fmt, county.fmt and user.fmt format library files. You can click on the '+' to the left of the format library file name to expand and view the formats contained in the file. The values and labels will be displayed on the right hand side of the screen when you select a format from the left. Only system administrators will be able to modify formats in the user.fmt format library.

To modify an existing format, click on the format on the left-hand side. You can add, delete or modify rows using the grid on the left-hand side. The buttons on the top of the screen update (or save) the format, insert a row, delete a row, or undo changes, respectively. The data type will determine what will be entered in the value field. The 'value' will be stored in the database. Labels will be used in all data entry screens, along with queries and reports. You can sort the entries by either value or label by right clicking on the column and selecting on "Sort by value" or "Sort by label". You should sort the format by label before saving the format. This will allow the system to display the formats properly for the 'type ahead' feature.

To add a new format, click on New Format. This will bring up a blank format. You will need to supply the format name, the data type, and all entries (value, label).

User-defined Fields

Introduction to User-defined Fields

User defined fields are a new feature in CaST. An administrator can set up fields to collect information outside of the MDE definition. Each table in the database has a corresponding user-defined table associated with it. User-defined fields will be displayed and available for data entry in the same location as the associated table fields. For example, user-defined fields in the Baseline_User table will be displayed as an additional tab in the Baseline Information data entry screen. The user-defined fields for the procedure tables will be at the end of the procedures grid data entry section. The maximum number of user-defined fields per table is controlled by Microsoft Access, and is currently 200 per table.

The system includes several user-defined fields included in the previous version of CaST. Administrators at existing CaST sites can rename these fields according to how they entered data in CaST 3.01. New users of CaST can delete these fields. These fields include:

User-defined Table	Field Name
BreastCycle_User	B_Com1
BreastProcedures_User	BS_L1
	BS_L2
	BS_L3
	BS_L4
CervicalCycle_User	C_Com1
CervicalProcedures_User	CS_L1
	CS_L2
	CS_L3
	CS_L4

Decentralized programs: administrators should setup all user-defined fields before the distribution of the software to the local sites.

User-defined Fields Setup

Important: Administrators should make a backup of the database and make sure that no one else is using the system before modifying user-defined fields.

To create user-defined fields, system administrators can select Program Defined Fields from the Administration menu. The following screen will

appear:

Field Name	Field Type	Field Length	Format
Patient_Id	Text	15	
B_Cycle_Number	Integer	4	
BS_Procedure_Number	Integer	4	
BS_L1	Integer	4	
BS_L2	Integer	4	
BS_L3	Text	15	
BS_L4	Text	10	
*			

To select the table you wish to modify, select the table name from the drop down labeled "User Defined Table". Once the table is selected, the current fields in the table are displayed. Key fields, such as Patient ID, cycle number and procedure number, are always in the tables, and can not be modified or deleted.

To insert or delete a row, use the buttons at the top of the screen. There is also a restructure and an Undo button at the top. The following describes the columns involved in creating a user-defined field.

Field_Name - name of the field for database - won't be displayed on data entry screens. Example: smoking_status. The field must start with a letter, can not contain spaces, and can not begin with the word 'Field'. **Note:** you can not use a field name that is already in use by another table. The system will warn you when you try to add a field name that already exists.

Field_Type - choose from Text, Integer, Number, Date/Time, Yes/No, or Memo (Note: memo fields are not recommended)

Field Length - field length, you can leave blank for date/time, yes/no, or memo

Format Name - enter name of format if applicable. If you'd like your field to have a 'drop down' list, specify a valid format name (either a system one or one you created using the format library)

Edit_Include - automatically set to 1: no need to change

Display Name - used to label the field. Note, this label should be fairly short. The intent is not to have a long question or description, only to provide an easier to read label.

Display Width - used to control width of field (number of pixels). Set to 100 and then adjust if necessary.

When finished, click on Restructure to save changes. Changes to the tables will now be available in the data entry screens and queries.

Note: the order in which the fields are saved will be used for display and data entry purposes. The 'layout' options for the procedure grids can be used to re-order/hide the fields for these screens, but the order specified here will determine the layout for the baseline and cycle tabs. You can use the mouse to re-order the fields in the restructure screen. Select the row containing the field you wish to move, hold the mouse button down, and move the field to the desired location.


Please also note, there are limitations on what you can restructure during one process. For instance, if you are renaming field names, you can not move or insert fields at the same time. You must first rename the field, restructure, and then insert or move the fields. The system will warn you if you try to do both at the same time.

Backup/Restore/Compact

Backing up the Database

You should backup your database frequently, especially if your DATA directory does not get backed up automatically through a local backup utility or a network function. Even if there are regularly scheduled backups on your hard drive or network, do not depend on these to backup your CaST data.

To backup the CaST data, click on Backup from the Tools menu. The following screen will appear.



You have the option of backing up both the administrative and patient data, along with the 'user.fmt' file, which stores all user-defined formats. By default, the backup will reside in the BACKUP directory under the 'cast files' directory specified during software installation. The file(s) will be named with the current date and the type of database being backed up. If you wish to backup to another location, you can specify the drive and directory.


Note: backing up to a different device (network hard drive, zip disk, etc.) is **highly** recommended.

Click on Backup to begin the process.

Restoring a Backup of the Database

If you perform regular backups of your data and encounter data loss or data corruption, users with 'administrator' permission can restore the database from a previously backed up version.

First, make sure no one else is using the system. Next, select Restore from the Tools menu. The following screen will appear.




You have the option of restoring the main database (patient data), the administrative database (containing query definitions and users), or the 'user formats'. Select the database to restore and hit Restore. If the backup is located in a location other than the default, select the appropriate drive letter to locate the files.

Compacting the Database

Make sure you have a good backup of the database before you compact.

Compacting the database will free up wasted space in the database and make the database run more efficiently. You should compact the database every month or so.

First, make sure no one else is using the system. Next, select Compact from the Tools menu. The following screen will appear.



Select Compact to run the compact utility.

Frequently Asked Questions (FAQ)

Known Issues and Comments from current Users

This section will highlight known issues in the system. The following issues were either discovered in our own testing or discovered during the installation and use of CaST. We welcome additional comments on the system, we only ask that you review this section to determine if your problem or comment has already been addressed.

Is there a billing component in CaST?

There is not a billing component in CaST. CaST is used by over 30 States, Territories and Tribes. Adding a billing component that would satisfy all Programs, with the limited funds available for development, is not feasible at this time.

Programs can create separate billing systems and link them with the CaST databases. The CaST system uses Microsoft Access databases for storing the patient data. There are already a few Programs that have successfully creating billing systems. Programs can link the databases for viewing only, or even write code to update or add information to the CaST databases.

If your Program is considering adding a billing system to CaST, we would highly recommend you contact the CaST Technical Support Team.

How can I contact other CaST Users?

If you have a question that you would like to ask other CaST users, you can use the CaSTHelp email account. Send your question to casthelp@nbccedp.org, and put the word 'ListServ' in the subject line. The message will be reviewed and forwarded to the Data Managers at all CaST Programs. Please include your contact information in your email, as Data Managers will reply directly to you.

What do 'Old CBE - Normal' and 'Old CBE - Abnormal' mean?

You will notice these two options for the 'Consultant-Repeat Clinical Breast Exam' procedure. The DOS based version of CaST only included 'Normal' and 'Abnormal' responses for this procedure. At the request of Programs, and to be consistent with the 'screening' CBE procedure results, the Model Clinical Categories were added as valid responses. These existing data could not be converted to the Model Clinical Categories, so the formats were renamed to 'Old CBE – Normal' and 'Old CBE – Abnormal'. This does not affect the MDE export, since diagnostic test results are not exported, only whether procedures were performed. You should collect the Model Clinical Categories on your data forms, the older formats will only be used as placeholders for older data.

What are the 'Old Location' fields?

The Cycle and Procedure data entry screens have 'Old Location' fields. The DOS based version of CaST collected location (cycle and screening) using free text fields. The Windows based version of CaST includes the Enrollment Site/Provider List Manager, which creates drop down lists and limits the entry of locations. Due to the many ways the free text fields can be entered, the existing data was not converted. Instead, the older data was converted to these 'old location' fields. After the conversion, Programs were expected to setup the provider list and enter information into the new 'drop down' fields. The 'old location' fields can not be modified. The 'old location' field in the procedure 'grid' can be hidden, please see "Layout files for the procedure data entry screens".

Why doesn't CaST support the use of numbers for data entry?

CaST, by default, supports the type-ahead feature found in most Windows applications. This widely used method eliminates the need for users to keep up with the numbers that are stored in the database, which get formatted for display and reporting.

If your program wishes to use numbers instead of the character based method, you will not be able to use the type-ahead feature to search for text. You will also have to make sure that your data entry forms and data entry staff is aware of all numbers in the system.

If you wish to use numbers instead of character formats, you will have to modify the cast.fmt file located in the DATA directory. You can modify this file using a text editor, or with a slight modification to the file, the format editor in CaST. You can not modify the value stored in the database, only the Label which is displayed on the data entry screens and queries. If you are interested in modifying this file, it is very important that you contact the CaST Technical Support Team at IMS.

Why 'CaST II' for the directory and logon screen?

The first two digits in the version number of CaST will always coincide with the version of the software. The second two digits represent the version of the MDEs that the system will export. For example, version 5.0.50 is version 5.0 of the software, with a MDE export of 5.0. "CaST II" is used to differentiate between the DOS and Windows versions. The directory structure will stay the same for future versions.

How can I print the MDE summary and error files?

The viewer in CaST currently does not have a print function. Future versions will allow you to print the summary and error files. These files can be printed using NotePad or a word processor. The files are stored in the MDE subfolder of the Cast support files directory. The default support files directory is "C:\Cast Files".

Why doesn't the cycle number database field match the record location in CaST?

This will affect you if you include the cycle number field (C_Cycle_Number or B_Cycle_Number) in a query, or if you are searching for records using the record ID field from the MDE edit listings. These internal cycle numbers may not match the location found when you are browsing cycle records, using the record locators at the bottom of the cycle data entry screen. The reason for this is that the system does not renumber cycle or procedure numbers when records get deleted. These numbers are internal to the system and are used to uniquely identify the record in the database. The internal cycle numbers, which will match the record ID fields, can be located on the top right on the cycle data entry screens.

How do I delete a procedure?

You can right click on a procedure to delete it. Select the procedure 'row' you wish to delete, use your right mouse button, and the select 'Delete'. Important: hitting the delete cycle button at the bottom of the screen will delete the entire cycle, including all procedures associated with the cycle.

How can I get an electronic list of my providers for use in other applications?

If you have Access 97 or another database/software package which can read Access 97 tables, the enrollment/provider list is contained in a table called "Provider_Sites".

If you do not have the ability to view this table, you can use the provider list export function. This function is used by programs utilizing decentralized data entry, but it can also be used to create a text file (fixed length), containing selected or all providers. This export function creates two text files, one which you specify the name (.txt) and a file describing the file layout (.zds). Both files can be viewed/edited using a text editor (Notepad, Wordpad, etc).

You can also utilize the advanced SQL system queries included in CaST. You can run these queries, and print or export the results.

Technical Specifications

Database Structures

CaST data, by default, are stored in Microsoft Access 97 tables. The user does have the option, though, to convert these databases to Access 2000.

CaST uses a relational database model, using Patient_ID, cycle and procedures number to link the tables. Please see the appendix of the printed documentation for a description of the database model and all database tables and fields. This appendix also includes documentation on the formats used to populate the data in the fields requiring formatted values.

Database Conversion

For programs already using CaST 4.1.1, please see Upgrading the databases for CaST 5.0.50. Detailed documentation on the conversion from version 3.01 to 4.0 is available upon request.

For programs converting from other systems, the printed documentation on database structures will provide you will instructions on how to convert data into CaST.


For questions regarding the conversion of your existing data, please contact the CaST Technical Support Team at IMS.

MDE Export Specifications

The MDE export program, for the most part, is a straight export to the MDE file definition. The export program does create separate records for Mams and Paps, and therefore sets all appropriate fields to unknown or "not-needed" when appropriate. For the most part, though, information in the system is not modified to eliminate errors/edits from the MDE edit program. For instance, if a field is not completed, a blank response will be reported in the MDEs, instead of Unknown. The only exception to this is when the value of "Pending" is valid. In these cases, blank results are coded to "Pending".

Documentation on the MDE export is available in the EXTRAS directory on the installation CD-ROM.

CaST 5.0 Data Structure



Baseline (Baseline/Contact Information Data Entry Tabs)					
Database Field Name	Description/Data Entry Field Label	Type/Width			Format Name/Specifications/Skip Patterns
<i>Patient_Id</i>	<i>Patient Identifier (auto generated)</i>	Character	15		see Appendix III - uses FIPS code, data entry site code and counter from TRACKER table
SSN	SS (Social Security) Number	Character	11		system automatically adds dashes in database field, 999-99-9999
Alt_Id	Alt (Alternate) ID	Character	15		used if woman doesn't have social security number (not in TAB order)
Alt_Id_Type	(Alternate ID) Type	Character	15		text to describe alternate id used (not in TAB order)
Name_Last	Last Name	Character	20		
Name_First	First (Name)	Character	15		
Name_Middle	Middle (Name)	Character	15		
Name_Maiden	Maiden (Name)	Character	20		
DateofBirth	Date of Birth	Date			
Race1	Race 1	Character	2		Race
Race2	Race 2	Character	2		Race
Race3	Race 3	Character	2		Race
Race4	Race 4	Character	2		Race
Race5	Race 5	Character	2		Race
Race6	Race 6	Character	2		Race
Ethnicity	Ethnicity	Character	1		Ethnicity
PlaceofBirth_City	POB (Place of Birth) City	Character	20		
PlaceofBirth_State	(Place of Birth) State (Postal Code)	Character	2		StateAbbrAbbr (can be found in cast.fmt)
PlaceofBirth_Country	(Place of Birth) Country	Character	20		autofilled with "USA" if one of 50 states
Enroll_Loc	Enrollment Location	Character	7		from Provider_Sites table, with Enrollment field set to yes
Enroll_Date	Enrollment Date	Date			
ChronicDisease_Number	CD (Chronic Disease) #	Character	9		
Chart_Number	Chart #	Character	9		
Comments	Comments about patient	Memo			
Primary_Street	Primary address - street	Character	45		
Primary_City	Primary address - city	Character	20		
Primary_State	Primary address - state (postal code)	Character	2		StateAbbrAbbr (can be found in cast.fmt)
Primary_Zip	Primary address - zip code	Character	10		
Primary_County	Primary address - County	Character	20		uses list of counties from county.fmt file (stored as text)
Phone_Day	Day telephone number	Character	20		include special characters if you wish
Phone_Nite	Night telephone number	Character	20		include special characters if you wish
Sec_Street	Secondary address - street	Character	45		
Sec_City	Secondary address - city	Character	20		
Sec_State	Secondary address - state (postal code)	Character	2		StateAbbrAbbr (can be found in cast.fmt)
Sec_Zip	Secondary address - zip code	Character	10		
PatientStatus	PatientStatus	Character	1		PatientStatus
DateofDeath	Date of Death	Date			complete only if Patient Status = deceased
<i>I_User_Id</i>	<i>Last user to update record</i>	Character	10		
<i>I_User_Date</i>	<i>Date last modified</i>	Date			
<i>I_Modify_Date</i>	<i>Date last modified (client)</i>	Date			
<i>I_Create_Date</i>	<i>Date record was created</i>	Date			
<i>Patient_Id_Old</i>	<i>Holds ID from prevision version of system</i>	Character	30		CaST 3.01 (DOS) and earlier
<i>Converted_Record</i>	<i>Flag for converted records</i>	Logical	1		set to 1 (yes) for CaST 3.01 data

Appendix II

CaST 5.0 Field Descriptions/Formats/Skip Patterns

MDE fields are in BOLD (and Red) - system/read-only/hidden fields are in Italics (and Blue)

Breast Cycle (Patient Eligibility and Clinical History, Cycle Disposition, and Breast Stage/TX Information)				
Database Field Name	Description/Data Entry Field Label	Type/Width		Format Name/Specifications/Skip Patterns
<i>Patient_ID</i>	<i>Patient Identifier (auto generated)</i>	Character	15	
<i>B_Cycle_Number</i>	<i>Cycle number (auto generated)</i>	Numeric		
B_Screen_Loc	Cycle location	Character	7	from Provider_Sites table, where Primary_Provider field set to yes
<i>B_Screen_Loc_Old</i>	<i>Cycle location (free text field from CaST 3.01)</i>	Character	15	from CaST 3.01 data - read only in system (use B_Screen_Loc instead)
B_Eligible_Income	Income eligible	Character	1	YesNoUnknown
B_Eligible_Medicare_Aid	Eligible for medicare/caid	Character	1	YesNoUnknown
B_Eligible_Insurance	Eligible for insurance	Character	1	YesNoUnknown
B_Suppress_Reminders	Suppress reminders for this client	Character	1	YesNo
B_Suppress_MDE	Suppress MDE record for this cycle	Character	1	YesNo
B_Prior_Mam	Prior Mammogram	Character	1	YesNoUnknown - if yes, complete Prior Mam documented and date
B_Prior_Mam_Documented	Prior Mam documented	Character	1	YesNoUnknown
B_Prior_Mam_Date_MDY	Prior Mammogram date	Date		use 06 or 15 for unknown month or day
B_Referral_Source	Referral Source	Character	2	ReferralSource
B_Symptoms	Client Symptoms	Character	1	YesNoUnknown
B_Workup_Planned	Planned Workup for Breast Cancer	Character	1	DXWU - if Yes, complete status of final diagnosis
B_Final_DX_Status	Status of Final Diagnosis	Character	1	FinalDiagStatus - if 'Complete', fill in Final Diagnosis
B_Final_Diagnosis	Final Diagnosis	Character	1	BrstFinalDiag
B_Closure_Date	Date of Final Diagnosis	Date		complete for all 'status of final diagnosis' results, except BLANK and 'pending'
B_Tumor_Stage	Tumor stage	Character	1	BrstTumorStage - complete if Final Diagnosis = 'Invasive Breast Cancer'
B_Tumor_Size	Tumor size (cm)	Numeric		use blank or 0 for unknown - complete if Final Diagnosis = 'Invasive Breast Cancer'
B_Tx_Status	Status of treatment	Character	1	TxStatus - must be completed for 'Invasive', 'DCIS' and 'CIS'. Can be completed for 'LCIS'
B_Tx_Date	Treatment date	Date		use date treatment was started, or administrative date for other non 'pending' results (lost, refused, etc)
<i>B_User_Id</i>	<i>Last user to update this BR record (auto generated)</i>	Character	10	
<i>B_Modified_Date</i>	<i>Date last modified (BR record)</i>	Date		
<i>B_Create_Date</i>	<i>Date BR record was created (auto generated)</i>	Date		
<i>B_Study_Mod_Date</i>	<i>Date last modified (any record in cycle) (auto generated)</i>	Date		
<i>B_Converted_Record</i>	<i>Flag for converted records (auto generated)</i>	Logical	1	set to 1 (yes) for CaST 3.01 data
Breast Procedure (Breast Cycle Procedures Performed 'Grid')				
Database Field Name	Description/Data Entry Field Label	Type/Width		Format Name/Specifications/Skip Patterns
<i>Patient_Id</i>	<i>Patient Identifier (auto generated)</i>	Character	15	
<i>B_Cycle_Number</i>	<i>Cycle number (auto generated)</i>	Numeric		
<i>BS_Procedure_Number</i>	<i>Procedure Number (within cycle number) (auto generated)</i>	Numeric		sequential number for each procedure in cycle
BS_Name	(Procedure) name	Character	2	BreastProcedureName
BS_Requested_Date	Date procedure was requested	Date		
BS_Loc	Procedure location	Character	7	from Provider_Sites table, where Provider field set to yes
<i>BS_Loc_Old</i>	<i>Location procedure performed (free text field from CaST 3.01)</i>	Character	15	from CaST 3.01 data - field is read only in system (use BS_Loc instead)
BS_Performed_Date	Date procedure was performed	Date		
BS_Paid	Was procedure paid for by NBCCEDP funds	Character	1	YesNoUnknown3
BS_Result	Result of the procedure	Character	2	Depends on BS_Name, see Result Formats
BS_Result_Date	Date of results	Date		
BS_CPTCode	CPT Code for billing purposes	Character	10	CPTCodesBrstFmt
BS_FUP_Recommended	Recommendations	Character	2	BreastProcedureRecommendations
BS_Short_FUP_Months	Months for short term follow up	Numeric		
BS_Results_Notify_Date	Date client is notified about procedure results	Date		
<i>BS_User_Id</i>	<i>Last user to update this record (auto generated)</i>	Character	10	
<i>BS_Modified_Date</i>	<i>Date last modified (BS record) (auto generated)</i>	Date		
<i>BS_Create_Date</i>	<i>Date record was created (auto generated)</i>	Date		
<i>BS_Converted_Record</i>	<i>Flag for converted records (auto generated)</i>	Logical	1	set to 1 (yes) for CaST 3.01 data

Appendix II

CaST 5.0 Field Descriptions/Formats/Skip Patterns

MDE fields are in BOLD (and Red) - system/read-only/hidden fields are in Italics (and Blue)

Cervical Cycle (Patient Eligibility and Clinical History, Cycle Disposition, and Cervical Stage/TX Information)				
Database Field Name	Description/Data Entry Field Label	Type/Width		Format Name/Specifications/Skip Patterns
<i>Patient_Id</i>	<i>Patient Identifier (auto generated)</i>	Character	15	
<i>C_Cycle_Number</i>	<i>Cycle number (auto generated)</i>	Numeric		
C_Screen_Loc	Cycle location	Character	7	from Provider_Sites table, where Primary_Provider field set to yes
<i>C_Screen_Loc_Old</i>	<i>Cycle location (free text field from CaST 3.01)</i>	Character	15	from CaST 3.01 data - read only in system (use C_Screen_Loc instead)
C_Eligible_Income	Income Eligible	Character	1	YesNoUnknown
C_Eligible_Medicare_Aid	Eligible for medicare/caid	Character	1	YesNoUnknown
C_Eligible_Insurance	Eligible for insurance	Character	1	YesNoUnknown
C_Suppress_Reminders	Suppress reminders for this client	Character	1	YesNo
C_Suppress_MDE	Suppress MDE record for this cycle	Character	1	YesNo
C_Prior_Pap	Prior Pap Test	Character	1	YesNoUnknown - if yes, complete Prior Pap documented and date
C_Prior_Pap_Documented	Prior Pap documented	Character	1	YesNoUnknown
C_Prior_Pap_Date_MDY	Prior Pap Test Date	Date		use 06 or 15 for unknown month or day
C_Referral_Source	Referral source	Character	2	ReferralSource
C_Workup_Planned	Planned Workup for Cervical Cancer	Character	1	DXWU - if Yes, complete status of final diagnosis
C_Final_DX_Status	Status of Final Diagnosis	Character	1	FinalDiagStatus - if 'Complete', fill in Final Diagnosis
C_Final_Diagnosis	Final Diagnosis	Character	1	CervFinalDiag
C_Final_Diagnosis_Other_Text	Text for Other Final Diagnosis	Character	50	complete only if Final diagnosis is 'other'
C_Closure_Date	Date of final diagnosis	Date		complete for all 'status of final diagnosis' results, except BLANK and 'pending'
C_Tumor_Stage	Stage	Character	1	CervTumorStage complete if Final Diagnosis = 'Invasive Cancer'
C_Tx_Status	Status of treatment	Character	1	TxStatus must be completed for 'Invasive', 'CIN II', 'CIN III' or High Grade'
C_Tx_Date	Treatment date	Date		use date treatment was started, or administrative date for other non 'pending' results (lost, refused, etc)
<i>C_User_Id</i>	<i>Last user to update this CR record (auto generated)</i>	Character	10	
<i>C_Modified_Date</i>	<i>Date last modified (CR record) (auto generated)</i>	Date		
<i>C_Create_Date</i>	<i>Date CR record was created (auto generated)</i>	Date		
<i>C_Study_Mod_Date</i>	<i>Date last modified (Any record in cycle) (auto generated)</i>	Date		
<i>C_Converted_Record</i>	<i>Flag for converted records (auto generated)</i>	Logical	1	set to 1 (yes) for CaST 3.01 data
Cervical Procedure (Cervical Cycle Procedures Performed 'Grid')				
Database Field Name	Description/Data Entry Field Label	Type/Width		Format Name/Specifications/Skip Patterns
<i>Patient_Id</i>	<i>Patient identifier (auto generated)</i>	Character	15	
<i>C_Cycle_Number</i>	<i>Cycle Number (auto generated)</i>	Numeric		
<i>CS_Procedure_Number</i>	<i>Procedure Number (within cycle number) (auto generated)</i>	Numeric		sequential number for each procedure in cycle
CS_Name	Procedure name	Character	2	CervicalStudyName
CS_Requested_Date	Date study was requested	Date		
CS_Loc	Procedure location	Character	7	from Provider_Sites table, where Provider field set to yes
<i>CS_Loc_Old</i>	<i>Location procedure performed (free text field from CaST 3.01)</i>	Character	15	from CaST 3.01 data - read only in system (use CS_Loc instead)
CS_Performed_Date	Date study was performed	Date		
CS_Adequacy	Adequacy of the specimen (if cs_name = 1 or 13)	Character	1	CervicalPAPadequacy
CS_SpecimenType	Specimen Type for 2001 Pap Test (cs_name = 13)	Character	1	CervicalSpecimenType
CS_Paid	Was procedure paid for by NBCCEDP funds	Character	1	YesNoUnknown3
CS_Result	Result of the procedure	Character	2	Depends on CS_Name, see Result Formats
CS_Other_Result_Text	Other results text	Character	50	free text for CS_Name = 1 or 13, CS_Result =7 (NO skip pattern)
CS_Result_Date	Date of results	Date		
CS_CPTCode	CPT Code for billing purposes	Character	10	CPTCodesCervFmt
CS_FUP_Recommended	Recommendations	Character	2	CervicalFUP
CS_Short_FUP_Months	Months for short term follow up	Numeric		
CS_Results_Notify_Date	Date client is notified about procedure results	Date		
<i>CS_User_Id</i>	<i>Last user to update this record (auto generated)</i>	Character	10	
<i>CS_Modified_Date</i>	<i>Date last modified (BS record) (auto generated)</i>	Date		
<i>CS_Create_Date</i>	<i>Date record was created (auto generated)</i>	Date		
<i>CS_Converted_Record</i>	<i>Flag for converted records (auto generated)</i>	Logical	1	set to 1 (yes) for CaST 3.01 data

Appendix II

CaST 5.0 Field Descriptions/Formats/Skip Patterns

MDE fields are in BOLD (and Red) - system/read-only/hidden fields are in Italics (and Blue)

CaST Database Conversion/Format Specifications

February 2003

For programs converting data from other systems, an in-depth understanding of the CaST data structures is needed for a successful data conversion.

This document, along with Appendix II, will provide information not evident from the relational database model (Appendix I).

Relational Database Model

The relational database model shows the relationship and links between the tables in the castmain database. The relational links are as follows:

Baseline	-->	BreastCycle, CervicalCycle: linked by Patient_ID (one to many)
BreastCycle	-->	BreastProcedures, CervicalCycle -> CervicalProcedures: linked by Patient_ID and Cycle_Number
All tables	-->	User tables (one to one) – i.e. each record in Baseline must have a corresponding record in Baseline_User

Patient ID field

Patient_ID is created using the following three sections:

- A three digit representation of program FIPS code. The FIPS code (or program code assigned by the CDC) is padded out with zeros, so Alabama's FIPS code of "01" becomes "001".
- The three digit data entry site code. For centralized data entry sites, use "001"
- A nine digit sequential number, unique for each data entry site.

Example Patient_ID numbers are: 001001000000001, 001001000000002, etc.

The **Tracker** table in the main database keeps track of the program code, data entry site, and next available sequential number. After you convert the data records you will need to set variables in this table based upon the largest Patient_ID used in the conversion. The following fields will have to be completed in the Tracker table:

Program:	three digit FIPS code (zero pad)
Site:	data entry site (use 001 for centralized sites)
Pri_pgm:	two digit FIPS code
Pri_pgm_code:	two letter program abbreviation, i.e. AL = Alabama
Id_counter:	largest number used in the Patient_ID field PLUS 1. If you converted 10000 records, i.e. last ID = 001001000010000, enter 10001 in this field.

Provider Sites table

The provider_sites table contains key information on enrollment sites and providers. There are three types of sites: enrollment, primary provider, and provider. These correspond to the Enrollment location field in the Baseline table, the Cycle Location fields in the Cycle tables, and the location fields in the Procedures tables, respectively. If you want to populate formatted values for these fields, you will need to also populate the provider_sites table. If it is too difficult to convert these locations (old data is all text), you can populate the corresponding “_Old” fields, which are read only in the system. You can then setup these sites and start collecting these data, converting the older data as time permits.

The format for site_num in the Provider_sites table is three digit data entry site (001 for central sites) and a four digit sequential number (i.e. 0001, 0002). The following additional fields should be completed:

Site_Name:	to identify the site
Site_City:	required for MDE export purposes (especially for enrollment sites)
Site_State:	
County:	using formats in county.fmt file. Two digit fips and three digit county FIPS code
Enrollment:	Is site an enrollment site? If yes, site will appear in enrollment location in baseline form: yes/no (-1 = yes, 0 = no)
Primary_Provider:	Is site a primary provider. If yes, site will appear in cycle location fields in cycle forms: yes/no (-1 = yes, 0 = no)
Provider:	Is site a provider. If yes, site will appear in location fields in procedure forms: yes/no (-1 = yes, 0 = no)

Note: there are three additional flags in the Provider_sites table: status_enrollment, status_pri_provider, and status_provider. These flags are automatically set when modifications are made to the three Yes/No provider type fields in CaST (see above). These flags are used to determine if a site is ‘active’, ‘inactive’ or ‘never active’. When the ‘type’ is changed from yes to no, CaST searches the database for existing records. If found, the ‘status’ flag is set to ‘inactive’. The site will then appear as ‘inactive’ during data entry.

There is a function in CaST called ‘re-sync’ in the Provider List menu off of the main menu. If you populate the provider table outside of CaST and do NOT have code to set these additional status flags, you should run this function after any updates.

Formats

The following pages contain the formats used in CaST 5.0. Most of these formats are contained in the default format library (cast.fmt), the first few listed are contained in the user-defined library (user.fmt) and can be expanded by the system administrator.

These formats are referenced in Appendix II – CaST 5.0 Field Descriptions/Formats.

[Format=Ethnicity] (in user.fmt)

1=Hispanic origin
2=Non-Hispanic origin
9=Unknown

[Format=Race] (in user.fmt)

1=White
2=Black
3=Asian
4=Pacific Islander
5=American Indian
6=Eskimo
7=Aleutian Islander
8=Other
9=Unknown
51=Native Hawaiian
98=Am.Ind/Esk./Aleut Isl. (3.0 data only)
97=Asian/Pac.Is. (3.0 data only)

[Format=ReferralSource] (in user.fmt)

1=Self
2=Provider
3=Other
4=Outreach
5=BCCP Reminder
9=Unknown

[Format=PatientStatus]

1=Active
2=Deceased
3=Inactive
4=Out of area
5=Temporarily Inactive

[Format=DXWU]

1=Planned
2=Not planned
3=Not yet determined

[Format=FinalDiagStatus]

1=Complete
2=Pending
3=Lost to follow-up
4=Refused
5=Deceased

[Format=TxStatus]

0=Refused by client
1=Lost to follow-up
2=Transportation problems
3=Financial problems
4=Not indicated/Not needed
5=Other problems
6=Client deceased
7=Treatment Started
8=Pending/Unknown

[Format=YesNo]

1=Yes

2=No

[Format=YesNoUnknown]

1=Yes

2=No

9=Unknown

[Format=YesNoUnknown3]

1=Yes

2=No

3=Unknown

[Format=BreastProcedureName]

1=Mammogram (initial)

2=Additional mam views

3=Clinical Breast Exam (CBE)

4=Ultrasound

5=Surgical consultation

6=Fine needle aspirate (FNA)

7=Biopsy

8=Consultant-Repeat Clinical Breast Exam

[Format=BreastMAMResult] (for BS_Name = 1, 2)

0=Refused

1=Negative

10=Result Pending

11=Result unk, presumed abn, non-program funded

12=Done recently elsewhere, non-funded

2=Benign finding

3=Probably benign

4=Susp abn (consider Bx)

5=Highly sugg of malig

6=Assessment is incomplete

7=Unsatisfactory

8=Not needed

9=Needed but not perf

98=Refused/Not Done (3.0 data only)

[Format=BreastCBEResult] (for BS_Name = 3)

0=Refused

2=Benign finding

30=Normal exam

31=Discrete palp mass - Susp for Ca

32=Discrete palp mass (Dx Benign)

33=Bloody/serous nipple disch

34=Nipple/areolar scaliness

35=Skin dimpling/retraction

36=Not done - Normal CBE in past 12 months

37=Not done - oth/unk reason

[Format=BreastUltraResult] (for BS_Name = 4)

0=Refused
15=Normal/No abnormality noted
16=Cystic mass
17=Susp for malign
18=Other benign abnormality
37=Not done - oth/unk reason
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=BreastConsultResult] (for BS_Name = 5)

0=Refused
19=No intervention at this time - routine FU
20=Short term FU
21=Biopsy/FNA recommended
37=Not done - oth/unk reason
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=BreastAspResult] (for BS_Name = 6)

0=Refused
22=No fluid/tissue obtained
23=Not Susp for Ca
24=Susp for ca
37=Not done - oth/unk reason
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=BreastBiopsyResult] (for BS_Name = 7)

0=Refused
25=Hyperplasia
26=Other benign chg
27=CIS - other
28=Invasive breast ca
29=Normal breast tissue
37=Not done - oth/unk reason
38=Ductal CIS
39=Lobular CIS
40=Atypical ductal hyperplasia (ADH)
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=BreastConsultCBEResult] (for BS_Name = 8)

0=Refused
13=Old CBE - Normal
14=Old CBE - Abnormal
2=Benign finding
30=Normal exam
31=Discrete palp mass - Susp for Ca
32=Discrete palp mass (Dx Benign)
33=Bloody/serous nipple disch
34=Nipple/areolar scaliness
35=Skin dimpling/retraction
37=Not done - oth/unk reason
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=BreastProcedureRecommendations]

- 1=Follow routine screening
- 10=Obtain definitive Rx
- 2=Short-term FU mam
- 3=Additional mam views
- 4=Repeat mam immed
- 5=CBE by consult
- 6=Ultrasound
- 7=Surgical consult
- 8=Fine needle aspirate (FNA)
- 9=Biopsy

[Format=BrstFinalDiag]

- 1=Carcinma In Situ - other (CIS) - Stage 0
- 2=Invasive Breast Cancer
- 3=Breast Cancer not diagnosed
- 4=Ductal Carcinma In Situ (DCIS) - Stage 0
- 5=Lobular Carcinoma In Situ (LCIS) - Stage 0

[Format=BrstTumorStage]

- 1=AJCC Stage I
- 2=AJCC Stage II
- 3=AJCC Stage III
- 4=AJCC Stage IV
- 5=Summary Local
- 6=Summary Regional
- 7=Summary Distant
- 9=Unknown
- 8=Unstaged

[Format=CervicalStudyName]

- 1=Pap test 1991
- 10=Gynecologic consultation
- 11=Pelvic Exam
- 12=Colposcopy without biopsy
- 13=Pap test
- 14=HPV test
- 2=Colposcopy with biopsy
- 3=Pelvic Ultrasound
- 4=Other biopsy - not colposcopic
- 5=Endocervical curettage
- 6=LEEP
- 7=Cone
- 8=Hysterectomy

[Format=CervicalPapResults] (for CS_Name = 1)

0=Refused
1=Negative (WNL)
10=Needed but not perf
11=Result pending
12=Result unk, presumed abn, non-program funded
13=Done recently elsewhere, non-funded
2=Inf/Inflam/React Changes
27=Not done - oth/unk reason
3=ASCUS
30=AGUS
31=Adenocarcinoma
4=Low grade SIL/HPV
5=High grade SIL
6=Squamous cell carcinoma
7=Other
8=Unsatisfactory
9=Not needed
98=Refused/Not Done (3.0 data only)

[Format=CervicalPapResults2001] (for CS_Name = 13)

31=Adenocarcinoma
30=AGUS
33=ASC-H
32=ASC-US
13=Done recently elsewhere, non-funded
5=High grade SIL
4=Low grade SIL/HPV
10=Needed but not perf
34=Negative for intra. lesion or malig.
27=Not done - oth/unk reason
9=Not needed
7=Other
0=Refused
11=Result pending
12=Result unk, presumed abn, non-program funded
6=Squamous cell carcinoma

[Format=CervicalSTDResults] (for CS_Name = 2, 4, 6, 7, 8)

31=Adenocarcinoma
14=CIN I
15=CIN II
17=CIN III/CIS
18=Invasive Carcinoma
1=Negative (WNL)
27=Not done - oth/unk reason
19=Oth nonmalignant abn (HPV,condyloma)
0=Refused
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=CervicalUltraResult] (for CS_Name = 3)

0=Refused
1=Negative (WNL)
20=Abn - not susp for ca
21=Abn - susp for ca
27=Not done - oth/unk reason
28=Abnormal Pelvic
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=CervicalEndoResults] (for CS_Name = 5)

0=Refused
1=Negative (WNL)
14=CIN I
15=CIN II
17=CIN III/CIS
18=Invasive Carcinoma
19=Oth nonmalignant abn (HPV,condyloma)
24=No cancer tissue
27=Not done - oth/unk reason
31=Adenocarcinoma
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=CervicalConsultResults] (for CS_Name = 10, 12)

0=Refused
1=Negative (WNL)
2=Inf/Inflam/React Changes
23=Other abnormality
27=Not done - oth/unk reason
8=Unsatisfactory
98=Refused/Not Done (3.0 data only)
99=Unknown

[Format=CervicalPelvicResults] (for CS_Name = 11)

0=Refused
20=Abn - not susp for ca
21=Abn - susp for ca
22=Not done - Normal PE in past 12 months
26=Normal
27=Not done - oth/unk reason
28=Abnormal Pelvic
29=Not indicated/Not needed
98=Refused/Not Done (3.0 data only)

[Format=CervicalHPVTestResults] (for CS_Name = 14)

35=Negative (HPV)
36=Positive (HPV)

[Format=CervicalPAPAdequacy] (for CS_Name = 1)

1=Satisfactory
2=Satisfactory for Eval, but limited by...
3=Unsatisfactory
9=Unknown

[Format=CervicalPAPAdequacy2001] (for CS_Name = 13)

1=Satisfactory
9=Unknown
3=Unsatisfactory

[Format=CervicalSpecimenType] (for CS_Name = 13)

1=Conventional smear
2=Liquid Based
3=Other
4=Unknown

[Format=CervFinalDiag]

1=Normal/Benign reaction/inflammation
2=HPV/Condylomata/Atypia
3=CINI/mild dysplasia (bx dx)
4=CINII/moderate dysplasia (bx dx)
5=CINIII/severe dysplasia/CIS (Stage 0) (bx dx)
6=Invasive Cervical Carcinoma (bx dx)
7=Other
8=Low Grade SIL (bx dx)
9=High Grade SIL (bx dx)

[Format=CervTumorStage]

1=Stage I
2=Stage II
3=Stage III
4=Stage IV
5=Summary Local
6=Summary Regional
7=Summary Distant
8=Unstaged
9=Unknown

[Format=CervicalFUP]

1=Follow routine screening
10=LEEP
11=Cone
12=HPV Test
2=Short-term follow up
3=Repeat Pap smear immediately
4=Colposcopy
5=Pelvic Ultrasound
6=Other biopsy
7=Gynecologic consultation
8=Definitive treatment
9=Hysterectomy

[Format=CPTCodesCervFmt] (in user.fmt)

57452=57452-Colposcopy without Biopsy

57454=57454-Colposcopy with Biopsy

88150=88150-Pap Smear

88151=88151-Repeat Pap Smear

88156=88156-Pap Smear, reported in Bethesda System

88157=88157-Pap smear, reported in the Bethesda System requiring interpretation by physician

88305=88305-Colposcopy Biopsy Interpretation

99201=99201-New Patient - Office Visit (10 minutes face to face)

99202=99202-New Patient - Office Visit (20 minutes face to face)

99203=99203-New Patient - Office Visit (30 minutes face to face)

99211=99211-Established Patient - Office Visit (05 minutes face to face)

99212=99212-Established Patient - Office Visit (10 minutes face to face)

99213=99213-Established Patient - Office Visit (15 minutes face to face)

99241=99241-Consultation Visit (15 minutes face to face)

99242=99242-Consultation Visit (30 minutes face to face)

99243=99243-Consultation Visit (40 minutes face to face)

[Format=CPTCodesBrstFmt] (in user.fmt)

19000=19000-Aspiration of Cyst of Breast

19001=19001-Aspiration of cyst of breast, additional

19100=19100-Biopsy of breast, needle core (surgical procedure only)

19101=19101-Incisional biopsy of breast

19120=19120-Excision of cyst, fibroadenoma, or other benign or malignant tumor aberrant breast tissue, duct lesion or nipple lesion

19125=19125-Excision of breast lesion identified by pre-operative placement of radiological marker - single lesion

19126=19126-Excision of breast lesion identified by pre-operative placement of radiological marker - each additional lesion

19290=19290-Preoperative placement of needle localization wire - breast

76090=76090-Diagnostic/Follow-up - Unilateral Mammogram

76091=76091-Diagnostic/Follow-up - Bilateral Mammogram

76092=76092-Screening mammogram

76095=76095-Stereotactic localization for breast biopsy, each lesion, radiological supervision and interpretation

76096=76096-Preoperative placement of needle localization wire, breast, radiological supervision and interpretation

76098=76098-Radiological examination, surgical specimen

76938=76938-Ultrasonic guidance for cyst aspiration, radiological supervision and interpretation

76942=76942-Ultrasonic guidance for needle biopsy, radiological supervision and interpretation

76645=76645-Ultrasound - Echography, Breasts

88170=88170-Fine Needle Aspiration with/without preparation of smears

88172=88172-Evaluation of Fine Needle Aspiration

88173=88173-Interpretation and Report of Fine Needle Aspiration

88305=88305-Breast biopsy interpretation

99201=99201-New Patient - Office Visit (10 minutes face to face)

99202=99202-New Patient - Office Visit (20 minutes face to face)

99203=99203-New Patient - Office Visit (30 minutes face to face)

99211=99211-Established Patient - Office Visit (05 minutes face to face)

99212=99212-Established Patient - Office Visit (10 minutes face to face)

99213=99213-Established Patient - Office Visit (15 minutes face to face)

99241=99241-Consultation Visit (15 minutes face to face)

99242=99242-Consultation Visit (30 minutes face to face)

99243=99243-Consultation Visit (40 minutes face to face)